

## PRODUCT INFORMATION

**Model: RAFR14041**

**Type: Labcold Advanced Laboratory Refrigerator**

For laboratory and general use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.



**Capacity: 385 litres (gross)**

**Construction: White enamel exterior, stainless steel interior. Cabinet mounted on castors.**

**Shelves supplied: 4; Plastic coated wire mesh**

**Exterior dimensions (mm): 1995 high x 600 wide x 700 deep**

### Features:

**Microprocessor temperature control with digital display**

**Audible and visible hi/lo temperature alarms**

**Fan circulated interior**

**Door lock**

**CFC and HCFC free**

**Temperature set-point range: +2 to +12°C**

**Automatic defrost with condensate removal**

**Maximum ambient temperature: 43°C**

**Order reference: RAFR 14041**

**N.B: Remote alarm contacts no longer supplied as standard but available at additional cost – order ref. REMCONADV**