

PRODUCT INFORMATION

Model: RLAA5002 (Factory fitted)

RLAA5002/1 (Stand-alone)

**Type: Labcold Digital Temperature Display/Alarm
With max/min record**

For use with any refrigerator or freezer operating in the range: +70°C to -50°C

Available either factory fitted in a specially constructed steel bracket (part ref RLAA5002); or as a stand-alone item (part ref RLAA5002/1)



RLAA 50002



As fitted to a Labcold Refrigerator

Temperature range: + 70°C to -50°C

Resolution: 1°C

Battery: 1 x AAA (supplied)

Exterior dimensions excluding bracket (mm): 115 high x 62 wide x 15 deep

Door or panel mounted as appropriate to cabinet construction.

Mounted in purpose built steel bracket (factory fitted version only)

Features:

Digital display of temperature

Audible alarm with mute facility

Maximum and minimum temperature displays with reset

Switchable between Fahrenheit and Centigrade

Low battery indicator

Detachable probe with 3.3m lead

Lead coiled within refrigeration chamber to allow user preference of monitoring position

Order reference: RLAA5002; RLAA5002/1

Labcold Limited, The Loddon Centre, Wade Road, Basingstoke, RG24 8FL, UK