

## **PRODUCT INFORMATION**

**Model: RLCH0100**

**Type: Labcold Cooled Incubator (Glass door)**

Benchtop controlled temperature incubator for maintaining samples at temperatures accurately controlled within the range +5 to +40°C.



**Capacity: 88 litres (gross)**

**Construction: White enamel exterior, Stainless steel interior, Double glazed door.**

**Shelves supplied: 2; Plastic coated wire mesh**

**Exterior dimensions (mm): 810H x 480W x 530D**

**Features:**

**Microprocessor temperature control with digital display**

**Temperature set-point range: +5 to +40°C**

**Temperature precision: better than  $\pm 0.5^\circ\text{C}$  at +5°C**

**Fan circulated interior**

**CFC and HFC free**

**Manual defrost with condensate drain**

**Maximum ambient temperature: +32°C**

**Order Reference: RLCH 0100**