

PRODUCT INFORMATION

Model: RLCH01501/AOP

Type: Labcold Cooled Incubator with Advanced Option Pack

For use in laboratories, hospitals, and other establishments to maintain samples, cultures and materials at a constant temperature. Suitable for bench use and can be used in conjunction with appropriate agitator (not included).



Capacity: 145 litres (gross)

Construction: White enamel exterior, Moulded plastic interior, Double glazed door.

Shelves supplied: 1; rigid steel stand to support agitator - removable

Exterior dimensions (mm): 1000H x 600W x 650D

Features:

Glass door and internal light

Fan circulated interior

Temperature range: +5 to +40°C (or can be factory preset)

Digital temperature controller and display

Temperature precision: better than $\pm 0.5^{\circ}\text{C}$ at +5°C

Labcold Alarm and Control System status indicator and alarm:

Audible and visible hi/lo temperature alarms

Door ajar and mains failure alarm

Remote alarm contacts

Seven day circular temperature chart recorder

Internal power socket

CFC and HFC free

Manual Defrost (required periodically) to maintain fixed temperature

Order Reference: RLCH 01501/AOP