

PRODUCT INFORMATION

Model: RLDF0105 Bench top/Wall mounted Refrigerator

Type: Labcold Pharmacy/Vaccine Refrigerator

A purpose-built refrigerator for laboratory, surgery, pharmacy and ward use for the storage of drugs and vaccines. Conforms fully with Department of Health “Green Book” requirements, plus recommendations in Royal Pharmaceutical Society of Great Britain’s “Code of Ethics and Standards”.



Capacity: 36 litres

Construction: White enamel exterior, metal interior

Shelves Supplied: 2 wire mesh

Exterior Dimensions (mm): 535H x 445W x 465D

Features:

Integrated Digital display of temperature with max/min record

Audible alarm with mute facility

Fan assisted air circulation

Integral wall mounting brackets

Door lock

CFC and HCFC free and Ammonia free

Accurate temperature control: +2°C to +8°C

Automatic defrost and removal of condensate

Order Reference: RLDF0105