

## PRODUCT INFORMATION

**Model: RLDG13043**

**Type: Labcold Pharmacy/Vaccine Refrigerator, Glass Door Model.**

For laboratory, surgery and general use for the storage of drugs and vaccines. Conforms in full with New UK Department of Health "Green Book" requirements.



**Capacity: 360 litres (gross)**

**Construction: White enamel exterior, fully moulded plastic interior.**

**Shelves Supplied: 6, Plastic coated wire mesh.**

**Exterior Dimensions (mm): 1825 high x 600 wide x 615 deep**

### Features:

- Double glazed door
- Digital display of temperature with max/min record
- Audible alarm with mute facility
- Fan assisted circulation
- Door Lock
- Easy-grip door handle
- CFC, HCFC and Ammonia free
- Accurate temperature control: +2°C to +8°C
- Automatic Defrost

**Order Reference: RLDG 13043**