

PRODUCT INFORMATION

Model: RLFF13246 and RLFF13246/LK (Lockable)

Type: Labcold Sparkfree Laboratory Fridge/Freezer

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



Capacities: Fridge section 230 Litres, Freezer section 160 Litres.
Construction: White enamel exterior, fully moulded plastic interior.
Internal configuration: Fridge: 4 Removable glass shelves.
Freezer: 4 Plastic drawers.
Exterior Dimensions (mm): 1965 high x 595 wide x 630 deep

Features:

- Spark free Interior (both compartments)**
- CFC and HCFC free**
- Temperature set-point range: +2 to +8°C (fridge); -17 to -23°C (freezer)**
- Automatic Defrost (fridge section only)**
- Removable door pockets**
- Independent temperature control of each compartment**

Please note: In order to provide a sparkfree interior no internal air circulating fan is fitted to this product. Some gradient in temperature, particularly in the refrigerator section, is to be expected.

Order Reference: RLFF 13246
RLFF 13426/LK