

PRODUCT INFORMATION

Model: RLPF09043

Type: Labcold Fan Assisted Laboratory Refrigerator

For laboratory and general use where improved temperature uniformity is required and there is no risk that volatile, flammable materials or components might be stored.



Image used for
illustration purposes

Capacity: 285 litres (gross)

Construction: White enamel exterior, fully moulded plastic interior

Shelves supplied: 5, glass shelves.

Exterior dimensions (mm): 1440 high x 600 wide x 635 deep

Features:

Fan assisted air circulation

Door lock

CFC and HCFC free

Temperature set-point range: +2 to +8°C

Automatic defrost

Order reference: RLPF 09043