

PRODUCT INFORMATION

Model: RLPL04043

Type: Labcold Basic Refrigerator

For laboratory and general use where reasonable performance and low cost are important and a sparkfree interior is not required.



Capacity: 105 Litres (gross)

Construction: White enamel exterior, moulded plastic interior

Shelves Supplied: 3 glass

Exterior Dimensions (mm): 845 high x 494 wide x 494 deep

Features:

Internal light and thermostat

CFC and HFC free

Temperature set-point range: +2 to +8°C

Automatic Defrost

Order Reference: RLPL 04043