

PRODUCT INFORMATION

Model: RLSD01501

Type: Labcold Refrigerated (Cooled) Incubator

For maintaining samples at temperatures accurately controlled within the range 0 to 40°C*.
Temperatures can be cycled where appropriate accessories are fitted.



Capacity: 145 litres (gross)

Construction: White enamel exterior, moulded plastic interior.

Shelves supplied: 3; Plastic coated wire mesh

Exterior dimensions (mm): 860 high x 600 wide x 650 deep

Features:

Microprocessor temperature control with digital display

Temperature set-point range: 0 to 40°C*

Temperature precision: better than +/- 0.5°C at 5°C

Fan circulated interior

CFC and HCFC free

Manual defrost with condensate drain

Maximum ambient temperature: 32°C

***Alternative temperature ranges available to special order**

Order reference: RLSD 01501