

## **PRODUCT INFORMATION**

**Model: RLVF1825**

**Type: Labcold Large Capacity Economy Freezer**

For laboratory and general use where large capacity combined with value for money is required and there is no risk that volatile, flammable materials or components might be stored. Also suitable for storage of certain temperature sensitive materials such as enzymes.



**Capacity: 505 litres (gross)**

**Construction: White enamelled steel exterior and interior.**

**Shelves supplied: 6; fixed position, refrigerated**

**Exterior dimensions (mm): 1710 high x 765 wide x 725 deep**

**Features:**

**Door lock**

**CFC and HCFC free**

**Temperature set-point range: -18 to -25°C**

**Analogue temperature display**

**High temperature warning light**

**Operator controlled defrost**

**Order reference: RLVF 1825**