

## **PRODUCT INFORMATION**

**Model: RLVL03203**

**Type: Labcold Basic Freezer**

For laboratory and general use where reasonable performance and low cost are important.



**Capacity: 70 Litres (gross)**

**Construction: White enamel exterior, fully moulded plastic interior.**

**Internal arrangement: Three drawers sections separated by fixed refrigerated shelves.**

**Exterior Dimensions (mm): 845 high x 494 wide x 494 deep**

### **Features:**

**Two removable freezer drawers**

**CFC and HCFC free**

**Temperature set-point range: -15 to -20°C**

**Operator controlled defrost**

**Order Reference: RLVL 03203**