

PRODUCT INFORMATION

Model: RPF44042

Type: Labcold Advanced Pharmacy Refrigerator (Glass door)

For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.



Capacity: 1245 litres (gross)

Construction: White and blue enamel exterior, stainless steel interior.

Cabinet mounted on castors, lockable.

Shelves Supplied: 8; Plastic coated wire mesh

Exterior Dimensions (mm): 1980 high x 1440 wide x 820 deep

Features:

Two double glazed doors

Microprocessor temperature control with digital display

Audible and visible hi/lo temperature alarms

Fan circulated interior

Door Lock

CFC and HCFC free

Approximate controllable range: +2 to +8°C

Accuracy: at +4°C ± 3°C, including defrost

Automatic defrost with condensate removal

Maximum ambient temperature: +32°C

Order Reference: RPF 44042

Optional Accessories: Remote Alarm Contacts: REMCONADV

32mm/50mm Access Ports: RLAP00*

Wall Mount Chart Recorder: RECW5050