

## PRODUCT INFORMATION

**Model: RPFR05042 (Solid door)**  
**RPFR05042/GD (Glass door)**  
**Type: Labcold Pharmacy Refrigerator**

For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.



**Capacity: 129 litres (gross)**  
**Construction: AISI 304 grade stainless steel interior and exterior.**  
**GD suffix: Denotes refrigerator is fitted with glass door.**  
**Shelves supplied: 3; Plastic coated wire mesh 530 x 325mm**  
**Exterior dimensions (mm): 850 high x 630 wide x 567 deep**

### Features:

- Microprocessor temperature control with digital display to 1°C
- Hi/Lo temperature alarm – audible
- Fan circulated interior
- CFC and HCFC free
- Temperature set-point range: +2 to +8°C
- Automatic defrost with condensate removal
- Maximum ambient temperature: +43°C
- Internal light on glass door models only

**Order reference: RPFR 0504**  
**RPFR 05042/GD**