

Customer information

09.06.2020, version 1, RaKr

Parameter update for compact photometer PF-12^{Plus}

Overview	
Photometer	Compact photometer PF-12 ^{Plus}
Parameter version	(9)
Availability	from 09 th June 2020
Price	free of charge
Link	https://www.mn-net.com/de/nanocolor-device-updates
Required accessories	USB-AB-Connection cable (REF 919687) Computer



Update procedure

- The parameter file of the PF-12^{Plus} is updated with the parameter update software. An instruction manual on using this software can be found in the download package, which is available via the above mentioned link.
- After a successful parameter update the following version number will be shown during the device start behind the firmware version in brackets: **(9)**

Modifications

- Implementation of the new VISOCOLOR[®] Powder Pillow test kit Iron (REF 936227)
- Bugfixes

Contact

If you have further questions, please do not hesitate to contact us:

Technical support and customer service	Tel.:	+49 24 21 969-332
	E-Mail:	csc@mn-net.com
For orders (international)	E-Mail:	sales@mn-net.com

www.mn-net.com



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: sales-us@mn-net.com