


# MS-H280-Pro Magnetic Hotplate Stirrer



The Camlab Choice MS-H280-Pro is a small bench top magnetic hotplate stirrer that heats up to 280°C. It has a 135mm stainless steel hotplate with excellent corrosion resistance. The unit can be fitted with an external temperature sensor PT1000 for real time accuracy of +/- 1°C. The temperature and speed are displayed clearly by the LED panel. HOT warning indicates residual hotplate temperature after use.


## Specifications


Part Number	1190949
Frequency	50/60Hz
Max stirring qty	3L
Speed range	100-1500 rpm
Heat output	500W
Safety temp	320°C
Dimensions	150x260x80
Weight	1.8Kg
Workplate dimension	135mm



## Features

- Stainless steel hotplate
- DC brushless motor
- Temp Range RT- 280°C
- Variable speed 100-1500rpm





Camlab Limited, Camlab House, Norman Way, Over, Cambridge, CB24 5WE  
Tel: 01954 233110 • Fax: 01954 233101  
email: [mailbox@camlab.co.uk](mailto:mailbox@camlab.co.uk) • [www.camlab.co.uk](http://www.camlab.co.uk)