


MS7-H550-Pro Magnetic Hotplate Stirrer



The Camlab Choice MS7-H550-Pro LCD digital hotplate stirrer has a 7" square work plate which is ideal for larger vessels and containers. Widely used in chemical synthesis, bio-pharmaceutical and clinical laboratories, the MS7-H550S has an LCD display which shows target and actual temperature. It has optional PT1000 temperature sensors which display the temperature of the solution being stirred. The glass ceramic coated hotplate demonstrates excellent corrosion resistance and is very easy to clean.


Specifications

Part Number	1191813
Frequency	50/60Hz
Max stirring qty	20L
Speed range	100-1500 rpm
Heat output	1000W
Safety temp	580°C
Dimensions	215x360x112
Weight	5.3Kg
Workplate dimension	184 x 184mm



Features

- Ceramic plate
- DC brushless motor
- Temp RT - 550°C
- Variable speed





camlab

Camlab Limited, Camlab House, Norman Way, Over, Cambridge, CB24 5WE
Tel: 01954 233110 • Fax: 01954 233101
email: mailbox@camlab.co.uk • www.camlab.co.uk