


MSE

Small Capacity Bench Top Centrifuges

- Simple - Easy operation
- Fast - Up to 14,000rpm
- Quiet - Low noise output
- Flexible - 3 acceleration and brake rates
- Practical - Just 5 kilogrammes
- Safe - Lid locking system
- Functional - Up to 24 microtubes per run
- Dependable - Reliable

The **Micro Centaur Plus** from MSE



MSE UK Ltd

Distributed by Camlab Ltd

Camlab Ltd, 24 Norman Way Industrial Estate,
Over, Cambridge CB24 5WE, UK

Tel: +44 1954 233 110
Fax: +44 1954 233 101
www.camlab.co.uk
sales@camlab.co.uk



Quality Laboratory Equipment

Micro Centaur Plus Technical Data

Height (mm)	180
Width (mm)	220
Depth (mm)	270
Weight (kg)	5
Maximum speed	14,000 rpm
Maximum RCF	14,243 g
Maximum capacity	24 x 2.2ml
Number of acceleration rates	3
Number of braking rates	3
Number of programmes	9
Pulse function	Yes
Timer	1-99 mins
Speed display	Digital LED \pm 10 rpm
Selectable temperature range	N/A
Catalogue number: 220/230V 50Hz	MSB010.CX3.5
Catalogue number: 120V 60Hz	MSB010.CX3.1

Tube or Sample Type	Rotor Part Number	Maximum Rotor Capacity
1.5ml tube	43117-615	24
2.2ml tube	43117-615	24
0.2ml PCR strips	43117-616	2 x 8 and 2 x 4
0.5ml tube	43117-617	18
50mm haematocrit capillaries	43117-618	24
0.4ml tube	43117-615*	24
0.5ml tube	43117-615*	12

* requires additional 81780.002 tube support

