

# Tube Rollers

Rolling Action for Complete Mixing



MX-T6-S



MX-T6-Pro

Camlab choice roller mixers are ideal for gentle mixing action which consists of a number of motor driven rollers which rotate at either a fixed or variable speed. When the samples are placed between the rollers in tubes or bottles typically, they are gently rolled. A rocking action is also applied whereby the rollers are raised and lowered at one end during rolling. This increases the effectiveness of the mixing.

## Applications:

Mixing of blood samples, viscous samples, and where minimum aeration is required or for aiding de-frosting of samples.

## Specifications

	<b>MX-T6-S</b>	<b>MX-T6-Pro</b>
Product Code	1191099	1191098
Speed range [rpm]	0-70	10-70
Power [W]	25	0
Max shaking weight [kg]	5	4
Speed display	Scale	LCD
Time setting range [min]	-	1-1199
Overall dimensions [LxWxHmm]	470x260x120	530x260x120
Permissible ambient temp [°C]	5-40	5-40
Table oscillation breaking device	Yes	Yes



## Features

- Rolling action for complete mixing
- Sturdy construction for high capacity
- Use with a wide variety of accessories
- MX-T6-Pro model features LCD display for speed and timer adjustment. Speed and timer can be viewed simultaneously.
- Electronic time switching clock for timing or continuous operation mode



camlab

[www.camlab.co.uk](http://www.camlab.co.uk)