

**NEW**  
ADDITIONAL MODULES  
FOR INCO

**Your preferred model!  
Choose it!**



**The new price-performance  
for CO<sub>2</sub> Incubators INCO:  
Lean perfection, top precision!**

UNBEATABLE  
**5 PROS**  
FOR THE NEW  
**INCO**

### 1 Excellent price-performance ratio

We never save on essentials! Even the price-attractive INCO standard version provides unique precision and reliability.

### 2 Fit for your application

As many functions as necessary, as much individuality as possible! Choose from six additional modules for your own special top model

**Comfort module:** two gas connections with quick release connectors, automatic switch-over of gas cylinders

**Hygiene module:** electropolished interior seamlessly welded by laser

**Communication module:** USB interface, „Celsius“ standard software for programming and logging, ring log memory, printer interface

**CO<sub>2</sub> module:** extended CO<sub>2</sub> range from 0 to 20 %

**Premium module:** includes Comfort, Hygiene, Communication and CO<sub>2</sub> module

**Humidity module:** active microprocessor-controlled humidification and dehumidification (40 - 97 % rh)

### 3 100% AtmoSAFE

Our INCO is 100 % AtmoSAFE! Absolutely precise control of temperature, humidity and CO<sub>2</sub> content guarantees a perfect atmosphere.

### 4 100% safe and hygienic

The working chamber including all fittings as well as all sensors is sterilisable at 160° C by an automatic 4-hours-programme.

### 5 Quick and flexible

Thanks to our flexible production and our modern logistic organisation also individual solutions can be produced quickly and flexible.

For more information and detailed technical data ask for our INCO brochure. We are looking forward to your contact under [sales@memmert.com](mailto:sales@memmert.com)!



[www.memmert.com](http://www.memmert.com)

Memmert GmbH + Co. KG  
P.O. Box 1720  
D-91107 Schwabach  
Tel. +49 (0) 9122 / 925 - 0  
Fax +49 (0) 9122 / 145 85  
E-Mail: [sales@memmert.com](mailto:sales@memmert.com)