

# Ohaus Centrifuge Overview

## Mini Centrifuges



Models #: **FCS306**  
**Speed Range**  
 0 rpm – 6,000 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 2,000 g  
**Maximum Capacity (Rotor)**  
 8 x 1.5 / 2.0 ml

## Micro Centrifuges



Models #: **FCS513 230V**  
**Speed Range**  
 200 rpm – 13,500 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 17,317 g  
**Maximum Capacity (Rotor)**  
 24 x 1.5 / 2.0 ml



Models #: **FCS515 230V**  
**Speed Range**  
 200 rpm – 15,200 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 21,953 g  
**Maximum Capacity (Rotor)**  
 44 x 1.5 / 2.0ml; 12 x 5 ml



Models #: **FCS513+R01 230V**  
**Speed Range**  
 200 rpm – 13,500 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 17,317 g  
**Maximum Capacity (Rotor)**  
 24 x 1.5 / 2.0 ml



Models #: **FCS515R 230V**  
**Speed Range**  
 200 rpm – 15,200 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 21,953 g  
**Maximum Capacity (Rotor)**  
 44 x 1.5 / 2.0ml; 12 x 5 ml

## Multi Purpose Centrifuges



Models #: **FCS707+R05 100-230V**  
**Speed Range**  
 200 rpm – 6,800 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 4,445 g  
**Maximum Capacity (Rotor)**  
 8 x 15 ml



Models #: **FCS706 230V**  
**Speed Range**  
 200 rpm – 6,000 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 4,427 g  
**Maximum Capacity (Rotor)**  
 6 x 50 ml

## Multi Purpose Pro Centrifuges



Models #: **FCS714 230V**  
**Speed Range**  
 200 rpm – 14,000 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 18,624 g  
**Maximum Capacity (Rotor)**  
 4 x 100 ml



Models #: **FCS718 230V**  
**Speed Range**  
 200 rpm – 18,000 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 23,542 g  
**Maximum Capacity (Rotor)**  
 4 x 100 ml



Models #: **FCS816 230V**  
**Speed Range**  
 200 rpm – 15,000 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 21,379 g  
**Maximum Capacity (Rotor)**  
 6 x 250 ml



Models #: **FCS718R 230V**  
**Speed Range**  
 200 rpm – 18,000 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 23,542 g  
**Maximum Capacity (Rotor)**  
 4 x 100 ml



Models #: **FCS816R 230V**  
**Speed Range**  
 200 rpm – 16,000 rpm  
**Maximum Relative Centrifuge Force (x g)**  
 24,325 g  
**Maximum Capacity (Rotor)**  
 6 x 250 ml