




Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207248	Description	Rotor, Angle, 24x1.5/2.0ml, ID, Sealable		
	Compatible with:	FC5714	FC5718	FC5718R	FC5816	FC5816R
	Camlab product Code	1206875	1206876	1206877	1206878	1206879
	max. speed	14000 rpm	14000 rpm	15000 rpm	15000 rpm	16000 rpm
	max. RCF	18624 xg	18624 xg	21379 xg	21379 xg	24325 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 24x1.5/2.0ml, ID, Sealable		1.5/2.0ml	11mm		24
1207351	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		24
1207352	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		24


Rotors and Accessories for Ohaus Frontier FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207267	Description	Rotor, Swing out, 4x250ml, ID		
	Compatible with:				FC5816	FC5816R
	Camlab product Code				1206878	1206879
	max. speed				4500 rpm	4500 rpm
	max. RCF				3780 xg	3780 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
1207338	Rack, 8x4-10ml D16mm Vac, 2/pk		4-10ml	16mm	Vac	32
1207366	Adapter, 9x1.5/2.0ml D11mm, 2/pk		1.5/2.0ml	11mm		36
1207373	Screw Cap, Rack of 250ml Rotor, 2/pk					
1207374	O-ring, Rack of 250ml Rotor, 10/pk					
1207323	Rack, 1x250ml D62mm FB, 2/pk		250ml	62mm	FB	4
1207324	Rack, 1x100ml D41mm RB, 2/pk		100ml	41mm	RB	4
1207325	Rack, 3x50ml D29mm RB, 2/pk		50ml	29mm	RB	12
1207326	Rack, 3x50ml D29mm FA, 2/pk		50ml	29mm	FA	12
1207327	Rack, 1x50ml D34mm RB, 2/pk, v3		50ml	34mm	RB	4
1207328	Rack, 2x50ml D29mm Rim, 2/pk		50ml	29mm	Rim	8
1207329	Rack, 1x50ml D29mm Rim, 2/pk, v3		50ml	29mm	Rim	4
1207330	Rack, 5x25ml D24.5mm RB, 2/pk		25ml	24.5mm	RB	20
1207331	Rack, 3x25ml D24.5mm Rim, 2/pk		25ml	24.5mm	Rim	12
1207332	Rack, 9x15ml D17mm RB, 2/pk		15ml	17mm	RB	36
1207333	Rack, 7x15ml D17mm FA, 2/pk		15ml	17mm	FA	28
1207334	Rack, 7x15ml D17mm Rim, 2/pk		15ml	17mm	Rim	28
1207335	Rack, 8x10ml D16mm RB, 2/pk		10ml	16mm	RB	32
1207336	Rack, 14x5-7ml D13mm RB, 2/pk		5-7ml	13mm	RB	56
1207337	Rack, 10x1.6-7ml D13mm Vac, 2/pk		1.6-7ml	13mm	Vac	40


Rotors and Accessories for Ohaus Frontier FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207249	Description	Rotor, Angle, 6x250ml FB, ID		
	Compatible with:				FC5816	FC5816R
	Camlab product Code				1206878	1206879
	max. speed				8000 rpm	8000 rpm
	max. RCF				10016 xg	10016 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No.of tubes/rotor
No Adapter	Rotor, Angle, 6x250ml FB, ID		250ml	62mm	FB	6
1207367	Adapter, 8x1.5ml D11mm, 2/pk		1.5ml	11mm		48
1207368	Adapter, 5x10ml D16mm RB, 2/pk		10ml	16mm	RB	30
1207369	Adapter, 4x15ml D17mm FA, 2/pk		15ml	17mm	FA	24
1207370	Adapter, 2x30ml D26mm RB, 2/pk		30ml	26mm	RB	12
1207371	Adapter, 1x50ml D29mm RB, 2/pk		50ml	29mm	RB	6
1207372	Adapter, 1x50ml D30mm FA, 2/pk		50ml	30mm	FA	6


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207270	Description	Rotor, Swing out, 2x3MTP, ID		
	Compatible with:	FC5714	FC5718	FC5718R	FC5816	FC5816R
	Camlab Product Code	1206875	1206876	1207877	1206878	1206879
	max. speed	4500 rpm	4500 rpm	4500 rpm	4500 rpm	4500 rpm
	max. RCF	2716 xg	2716 xg	2716 xg	2716 xg	2716 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Swing out, 2x3MTP, ID		MTP			6
1207288	Bucket, 3xMTP, 2/pk		MTP			6
1207356	Adapter, 48x1.5/2.0ml D11mm, 2/pk		1.5/2.0ml	11mm		96


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code 1207251	Description	Rotor, Angle, 6x85ml RB, ID			
	Compatible with:	FC5718	FC5718R	FC5816	FC5816R	
	Camlab product Code max. speed	1206876	1206877	1206878	1206879	
	max. RCF	9000 rpm	9000 rpm	11000 rpm	13000 rpm	
		10413 xg	10413 xg	15555 xg	21726 xg	
Accessories						
Camlab product Code	Description	Vol./Tube	Diameter	Tube type	No.of tubes/rotor	
No Adapter	Rotor, Angle, 6x85ml RB, ID	85ml	38mm	RB	6	
1207357	Adapter, 1x15ml D17mm RB, 2/pk	15ml	17mm	RB	6	
1207358	Adapter, 1x30ml D25mm RB, 2/pk	30ml	25mm	RB	6	
1207359	Adapter, 1x50ml D29mm RB, 2/pk, v2	50ml	29mm	RB	6	
1207360	Adapter, 1x50ml D29mm FA, 2/pk	50ml	29mm	FA	6	
1207361	Adapter, 1x15ml D17mm FA, 2/pk	15ml	17mm	FA	6	
1207362	Adapter, 2x15ml D17mm RB, 2/pk	15ml	17mm	RB	12	
1207363	Adapter, 1x16ml D18mm RB, 2/pk	16ml	18mm	RB	6	


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo	Camlab product Code	1207252	Description	Rotor, Angle, 4x85ml RB, ID, Hi		
	Compatible with:		FC5718	FC5718R	FC5816	FC5816R
	Camlab product Code		1206876	1206877	1206878	1206879
	max. speed		12000 rpm	12000 rpm	12000 rpm	12000 rpm
	max. RCF		14809 xg	14809 xg	14809 xg	14809 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 4x85ml RB, ID, Hi		85ml	38mm	RB	4
1207357	Adapter, 1x15ml D17mm RB, 2/pk		15ml	17mm	RB	4
1207358	Adapter, 1x30ml D25mm RB, 2/pk		30ml	25mm	RB	4
1207359	Adapter, 1x50ml D29mm RB, 2/pk, v2		50ml	29mm	RB	4
1207360	Adapter, 1x50ml D29mm FA, 2/pk		50ml	29mm	FA	4
1207361	Adapter, 1x15ml D17mm FA, 2/pk		15ml	17mm	FA	4
1207362	Adapter, 2x15ml D17mm RB, 2/pk		15ml	17mm	RB	8
1207363	Adapter, 1x16ml D18mm RB, 2/pk		16ml	18mm	RB	4


Rotors and Accessories for Ohaus Frontier FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207271	Description	Rotor, Swing out, 16x50ml, ID		
	Compatible with:				FC5816	FC5816R
	Camlab product Code				1206878	1206879
	max. speed				4500 rpm	4500 rpm
	max. RCF				3735 xg	3735 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
1207339	Rack, 4x50ml D29mm FA, 2/pk		50ml	29mm	FA	16
1207340	Rack, 10x15ml D17mm FA, 2/pk		15ml	17mm	FA	40
1207341	Rack, 14x5-7ml D13mm Vac/Sar, 2/pk		5-7ml	13mm	Vac/Sar	56
1207342	Rack, 12x4-10ml D16mm Vac/Sar, 2/pk		4-10ml	16mm	Vac/Sar	48


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

	Camlab product Code	1207253	Description	Rotor, Angle, 10x50ml FA, ID		
	Compatible with:		FC5718	FC5718R	FC5816	FC5816R
	Camlab product Code		1206876	1206877	1206878	1206879
	max. speed		7500 rpm	7500 rpm	9000 rpm	10500 rpm
	max. RCF		8174 xg	8174 xg	11771 xg	16022 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 10x50ml FA, ID		50ml	29mm	FA	10
1207353	Adapter, 1x1.5/2.0ml D11mm, 6/pk		1.5/2.0ml	11mm		10
1207345	Adapter, 1x7ml D13.5mm RB, 2/pk		7ml	13.5mm	RB	10
1207346	Adapter, 1x5ml D13.5mm, 2/pk		5ml	13.5mm		10
1207350	Adapter, 1x15ml D17mm RB, 2/pk, v2		15ml	17mm	RB	10


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207255	Description	Rotor, Angle, 6x50ml RB, ID, Hi		
	Compatible with		FC5718	FC5718R	FC5816	FC5816R
	Camlab product Code		1206876	1206877	1206878	1206879
	max. speed		12000 rpm	12000 rpm	13000 rpm	13000 rpm
			max. RCF			
			13522 xg	13522 xg	15869 xg	15869 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 6x50ml RB, ID, Hi		50ml	29mm	RB	6
1207347	Adapter, 1x30ml D26mm, 2/pk		30ml	26mm		6
1207348	Adapter, 1x16ml D18mm, 2/pk		16ml	18mm		6
1207349	Adapter, 1x15ml D17mm, 2/pk		15ml	17mm		6
1207364	Adapter, 1x15ml D17.5mm FA, 2/pk		15ml	17.5mm	FA	6


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207256	Description	Rotor, Angle, 30x15ml RB/FA, ID		
	Compatible with:	FC5714	FC5718	FC5718R	FC5816	FC5816R
	Camlab Product Code	1206875	1206876	1206877	1206878	1206879
	max. speed	4500 rpm	4500 rpm	4500 rpm	4500 rpm	4500 rpm
	max. RCF	2830 xg	2830 xg	2830 xg	2830 xg	2830 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No.of tubes/rotor
No Adapter	Rotor, Angle, 30x15ml RB/FA, ID		15ml	17mm	RB/FA	30
1207353	Adapter, 1x1.5/2.0ml D11mm, 6/pk		1.5/2.0ml	11mm		30
1207345	Adapter, 1x7ml D13.5mm RB, 2/pk		7ml	13.5mm	RB	30
1207346	Adapter, 1x5ml D13.5mm, 2/pk		5ml	13.5mm		30


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207257	Description	Rotor, Angle, 20x10ml RB, ID, Hi		
	Compatible with:		FC5718	FC5718R	FC5816	FC5816R
	Camlab Product Code		1206876	1206877	1206878	1206879
	max. speed		12000 rpm	12000 rpm	12000 rpm	12000 rpm
	max. RCF		15775 xg	15775 xg	15775 xg	15775 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 20x10ml RB, ID, Hi		10ml	16mm	RB	20

Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo	Camlab product Code	1207259	Description	Rotor, Angle, 44x1.5/2.0ml, ID		
	Compatible with:		FC5718	FC5718R	FC5816	FC5816R
	Camlab product Code		1206876	1206877	1206878	1206879
	max. speed		13500 rpm	13500 rpm	15000 rpm	16000 rpm
	max. RCF		17113 xg	17113 xg	21128 xg	24039 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 44x1.5/2.0ml, ID		1.5/2.0ml	11mm		44
1207351	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		44
1207352	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		44

Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207261	Description	Rotor, Angle, 24x1.5/2.0ml, ID		
	Compatible with:	FC5714	FC5718	FC5718R	FC5816	FC5816R
	Camlab Product Code	1206875	1206876	1206877	1206878	1206879
	max. speed	14000 rpm	14000 rpm	15000 rpm	15000 rpm	16000 rpm
	max. RCF	18624 xg	18624 xg	21379 xg	21379 xg	24325 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No.of tubes/rotor
No Adapter	Rotor, Angle, 24x1.5/2.0ml, ID		1.5/2.0ml	11mm		24
1207351	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		24
1207352	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		24