



Rotors and Accessories for Ohaus Frontier FC5714,FC5718 FC5718R, FC5816, FC5816R

	Camlab product Code	1207248	Description	Rotor, Angle, 24x1.5/2.0ml, ID, Sealable		
	Compatible with:	FC5714	FC5718	FC5718R	FC5816	FC5816R
	Camlab product Code	1206875	1206876	1206877	1206878	1206879
	max. speed	14000 rpm	14000 rpm	15000 rpm	15000 rpm	16000 rpm
	max. RCF	18624 xg	18624 xg	21379 xg	21379 xg	24325 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No.of tubes/rotor
No Adapter	Rotor, Angle, 24x1.5/2.0ml, ID, Sealable		1.5/2.0ml	11mm		24
1207351	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		24
1207352	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		24


Rotors and Accessories for Ohaus Frontier FC 5714, FC5718, FC5718R

	Camlab product Code	1207268	Description	Rotor, Swing out, 4x100ml, ID		
	Compatible with:	FC5714	FC5718	FC5718R		
	Camlab product Code	1206875	1206876	1206877		
	max. speed	4500 rpm	4500 rpm	4500 rpm		
	max. RCF	3350 xg	3350 xg	3350 xg		
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
1207289	Rack, 1x100ml D46mm RB, 2/pk		100ml	46mm	RB	4
1207290	Rack, 1x85ml D38mm RB, 2/pk		85ml	38mm	RB	4
1207291	Rack, 1x50ml D34mm RB, 2/pk		50ml	34mm	RB	4
1207292	Rack, 1x50ml D29mm RB, 2/pk		50ml	29mm	RB	4
1207293	Rack, 1x50ml D29mm FA, 2/pk		50ml	29mm	FA	4
1207294	Rack, 1x50ml D29mm Rim, 2/pk		50ml	29mm	Rim	4
1207295	Rack, 1x30ml D25mm Rim, 2/pk		30ml	25mm	Rim	4
1207296	Rack, 1x25ml D25mm RB, 2/pk		25ml	25mm	RB	4
1207297	Rack, 2x15ml D17mm FA, 2/pk		15ml	17mm	FA	8
1207298	Rack, 1x15ml D17mm Rim, 2/pk		15ml	17mm	Rim	4
1207299	Rack, 7x15ml D17mm RB, 2/pk		15ml	17mm	RB	28
1207300	Rack, 9x5-7ml D13mm RB, 2/pk		5-7ml	13mm	RB	36
1207301	Rack, 10x1.5ml D11mm, 2/pk		1.5ml	11mm		40
1207302	Rack, 9x4.5ml D12mm FB, 2/pk		4.5ml	12mm	FB	36
1207303	Rack, 7x5-7ml D13mm RB, 2/pk		5-7ml	13mm	RB	28
1207304	Rack, 7x5-7ml D13mm Vac, 2/pk		5-7ml	13mm	Vac	28
1207305	Rack, 5x4-10ml D16mm Sar, 2/pk		4-10ml	16mm	Sar	20


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

Rotor Photo 	Camlab product Code	1207269	Description	Rotor, Swing out, 4x100ml, ID, Sealable		
	Compatible with:	FC5714	FC5718	FC5718R		
	Camlab product Code	1206875	1206876	1207877		
	max. speed	4000 rpm	5000 rpm	5000 rpm		
	max. RCF	2486 xg	3885 xg	3885 xg		
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
1207283	Bucket, 100ml w/o Cap, 2/pk					
1207284	Bucket, 100ml w/ Cap, Sealable, 2/pk					
1207375	Screw Cap, 100ml Bucket, Sealable, 2/pk					
1207376	O-ring, 100ml Bucket, Sealable, 10/pk					
1207306	Rack, 1x100ml D40mm RB, 2/pk		100ml	40mm	RB	4
1207307	Rack, 1x85ml D38mm RB, 2/pk, v2		85ml	38mm	RB	4
1207308	Rack, 1x50ml D34mm RB, 2/pk, v2		50ml	34mm	RB	4
1207309	Rack, 1x50ml D29mm RB, 2/pk, v2		50ml	29mm	RB	4
1207310	Rack, 1x50ml D29mm FA, 2/pk, v2		50ml	29mm	FA	4
1207311	Rack, 1x50ml D29mm Rim, 2/pk, v2		50ml	29mm	Rim	4
1207312	Rack, 1x30ml D25mm Rim, 2/pk, v2		30ml	25mm	Rim	4
1207313	Rack, 1x25ml D25mm RB, 2/pk, v2		25ml	25mm	RB	4
1207314	Rack, 1x15ml D17mm FA, 2/pk		15ml	17mm	FA	4
1207315	Rack, 2x15ml D17mm FA, 2/pk, v2		15ml	17mm	FA	8
1207316	Rack, 2x15ml D17mm Rim, 2/pk		15ml	17mm	Rim	8
1207317	Rack, 4x15ml D17mm RB, 2/pk		15ml	17mm	RB	16
1207318	Rack, 7x5-7ml D13mm RB, 2/pk, v2		5-7ml	13mm	RB	28
1207319	Rack, 5x1.5ml D11mm, 2/pk		1.5ml	11mm		20
1207320	Rack, 3x10ml D17mm Sar, 2/pk		10ml	17mm	Sar	12
1207321	Rack, 4x5-7ml D13mm RB, 2/pk		5-7ml	13mm	RB	16
1207322	Rack, 4x10ml D17mm Vac, 2/pk		10ml	17mm	Vac	16
1207285	Bucket, 2x50ml D29mm FA, w/ Rack, 2/pk		50ml	29mm	FA	8
1207286	Bucket, 3x15ml D17mm FA, w/ Rack, 2/pk		15ml	17mm	FA	12
1207287	Cyto, Rectangular Bucket, 2/pk					
1207377	Cyto, Slide Carrier, 1/pk					
1207378	Cyto, Container, 2x1.5ml Cyl., 1/pk					
1207379	Cyto, Container, 2x3ml Conical, 1/pk					
1207380	Cyto, Container Cover, 1/pk					
1207381	Cyto, Filter Card, Cyl., 100/pk					
1207382	Cyto, Filter Card, Conical, 100/pk					


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207270	Description	Rotor, Swing out, 2x3MTP, ID		
	Compatible with:	FC5714	FC5718	FC5718R	FC5816	FC5816R
	Camlab Product Code	1206875	1206876	1207877	1206878	1206879
	max. speed	4500 rpm	4500 rpm	4500 rpm	4500 rpm	4500 rpm
	max. RCF	2716 xg	2716 xg	2716 xg	2716 xg	2716 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Swing out, 2x3MTP, ID		MTP			6
1207288	Bucket, 3xMTP, 2/pk		MTP			6
1207356	Adapter, 48x1.5/2.0ml D11mm, 2/pk		1.5/2.0ml	11mm		96


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R

Rotor Photo	Camlab product Code	1207250	Description	Rotor, Angle, 6x85ml RB, ID, Hi		
	Compatible with		FC5718	FC5718R		
	Camlab product Code		1206876	1206877		
	max. speed		11000 rpm	13500 rpm		
	max. RCF		13932 xg	20984 xg		
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 6x85ml RB, ID, Hi		85ml	38mm	RB	6
1207357	Adapter, 1x15ml D17mm RB, 2/pk		15ml	17mm	RB	6
1207358	Adapter, 1x30ml D25mm RB, 2/pk		30ml	25mm	RB	6
1207359	Adapter, 1x50ml D29mm RB, 2/pk, v2		50ml	29mm	RB	6
1207360	Adapter, 1x50ml D29mm FA, 2/pk		50ml	29mm	FA	6
1207361	Adapter, 1x15ml D17mm FA, 2/pk		15ml	17mm	FA	6
1207362	Adapter, 2x15ml D17mm RB, 2/pk		15ml	17mm	RB	12
1207363	Adapter, 1x16ml D18mm RB, 2/pk		16ml	18mm	RB	6


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code 1207251	Description	Rotor, Angle, 6x85ml RB, ID				
	Compatible with:	FC5718	FC5718R	FC5816	FC5816R		
	Camlab product Code max. speed	1206876	9000 rpm	9000 rpm	11000 rpm	13000 rpm	
	max. RCF		10413 xg	10413 xg	15555 xg	21726 xg	
Accessories							
Camlab product Code	Description	Vol./Tube	Diameter	Tube type	No.of tubes/rotor		
No Adapter	Rotor, Angle, 6x85ml RB, ID	85ml	38mm	RB	6		
1207357	Adapter, 1x15ml D17mm RB, 2/pk	15ml	17mm	RB	6		
1207358	Adapter, 1x30ml D25mm RB, 2/pk	30ml	25mm	RB	6		
1207359	Adapter, 1x50ml D29mm RB, 2/pk, v2	50ml	29mm	RB	6		
1207360	Adapter, 1x50ml D29mm FA, 2/pk	50ml	29mm	FA	6		
1207361	Adapter, 1x15ml D17mm FA, 2/pk	15ml	17mm	FA	6		
1207362	Adapter, 2x15ml D17mm RB, 2/pk	15ml	17mm	RB	12		
1207363	Adapter, 1x16ml D18mm RB, 2/pk	16ml	18mm	RB	6		


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo	Camlab product Code	1207252	Description	Rotor, Angle, 4x85ml RB, ID, Hi		
	Compatible with:		FC5718	FC5718R	FC5816	FC5816R
	Camlab product Code		1206876	1206877	1206878	1206879
	max. speed		12000 rpm	12000 rpm	12000 rpm	12000 rpm
	max. RCF		14809 xg	14809 xg	14809 xg	14809 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 4x85ml RB, ID, Hi		85ml	38mm	RB	4
1207357	Adapter, 1x15ml D17mm RB, 2/pk		15ml	17mm	RB	4
1207358	Adapter, 1x30ml D25mm RB, 2/pk		30ml	25mm	RB	4
1207359	Adapter, 1x50ml D29mm RB, 2/pk, v2		50ml	29mm	RB	4
1207360	Adapter, 1x50ml D29mm FA, 2/pk		50ml	29mm	FA	4
1207361	Adapter, 1x15ml D17mm FA, 2/pk		15ml	17mm	FA	4
1207362	Adapter, 2x15ml D17mm RB, 2/pk		15ml	17mm	RB	8
1207363	Adapter, 1x16ml D18mm RB, 2/pk		16ml	18mm	RB	4


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo	Camlab product Code	1207253	Description	Rotor, Angle, 10x50ml FA, ID		
	Compatible with:		FC5718	FC5718R	FC5816	FC5816R
	Camlab product Code		1206876	1206877	1206878	1206879
	max. speed		7500 rpm	7500 rpm	9000 rpm	10500 rpm
	max. RCF		8174 xg	8174 xg	11771 xg	16022 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 10x50ml FA, ID		50ml	29mm	FA	10
1207353	Adapter, 1x1.5/2.0ml D11mm, 6/pk		1.5/2.0ml	11mm		10
1207345	Adapter, 1x7ml D13.5mm RB, 2/pk		7ml	13.5mm	RB	10
1207346	Adapter, 1x5ml D13.5mm, 2/pk		5ml	13.5mm		10
1207350	Adapter, 1x15ml D17mm RB, 2/pk, v2		15ml	17mm	RB	10


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

Rotor Photo 	Camlab product Code	1207254	Description	Rotor, Angle, 6x50ml RB/FA, ID		
	Compatible with:	FC5714	FC5718	FC5718R		
	Camlab product Code	1206875	1206876	1206877		
	max. speed	6000 rpm	6000 rpm	6000 rpm		
	max. RCF	4427 xg	4427 xg	4427 xg		
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 6x50ml RB/FA, ID		50ml	29mm	RB/FA	6
1207353	Adapter, 1x1.5/2.0ml D11mm, 6/pk		1.5/2.0ml	11mm		6
1207345	Adapter, 1x7ml D13.5mm RB, 2/pk		7ml	13.5mm	RB	6
1207346	Adapter, 1x5ml D13.5mm, 2/pk		5ml	13.5mm		6
1207347	Adapter, 1x30ml D26mm, 2/pk		30ml	26mm		6
1207348	Adapter, 1x16ml D18mm, 2/pk		16ml	18mm		6
1207349	Adapter, 1x15ml D17mm, 2/pk		15ml	17mm		6
1207350	Adapter, 1x15ml D17mm RB, 2/pk, v2		15ml	17mm	RB	6


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207255	Description	Rotor, Angle, 6x50ml RB, ID, Hi		
	Compatible with		FC5718	FC5718R	FC5816	FC5816R
	Camlab product Code		1206876	1206877	1206878	1206879
	max. speed		12000 rpm	12000 rpm	13000 rpm	13000 rpm
	max. RCF		13522 xg	13522 xg	15869 xg	15869 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 6x50ml RB, ID, Hi		50ml	29mm	RB	6
1207347	Adapter, 1x30ml D26mm, 2/pk		30ml	26mm		6
1207348	Adapter, 1x16ml D18mm, 2/pk		16ml	18mm		6
1207349	Adapter, 1x15ml D17mm, 2/pk		15ml	17mm		6
1207364	Adapter, 1x15ml D17.5mm FA, 2/pk		15ml	17.5mm	FA	6


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207256	Description	Rotor, Angle, 30x15ml RB/FA, ID		
	Compatible with:	FC5714	FC5718	FC5718R	FC5816	FC5816R
	Camlab Product Code	1206875	1206876	1206877	1206878	1206879
	max. speed	4500 rpm	4500 rpm	4500 rpm	4500 rpm	4500 rpm
	max. RCF	2830 xg	2830 xg	2830 xg	2830 xg	2830 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No.of tubes/rotor
No Adapter	Rotor, Angle, 30x15ml RB/FA, ID		15ml	17mm	RB/FA	30
1207353	Adapter, 1x1.5/2.0ml D11mm, 6/pk		1.5/2.0ml	11mm		30
1207345	Adapter, 1x7ml D13.5mm RB, 2/pk		7ml	13.5mm	RB	30
1207346	Adapter, 1x5ml D13.5mm, 2/pk		5ml	13.5mm		30


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207257	Description	Rotor, Angle, 20x10ml RB, ID, Hi		
	Compatible with:		FC5718	FC5718R	FC5816	FC5816R
	Camlab Product Code		1206876	1206877	1206878	1206879
	max. speed		12000 rpm	12000 rpm	12000 rpm	12000 rpm
	max. RCF		15775 xg	15775 xg	15775 xg	15775 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 20x10ml RB, ID, Hi		10ml	16mm	RB	20


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

Rotor Photo 	Camlab product Code	1207258	Description	Rotor, Angle, 12x15ml RB/FA, ID		
	Compatible with:	FC5714	FC5718	FC5718R		
	Camlab Product Code	1206875	1206876	1206877		
	max. speed	6000 rpm	6000 rpm	6000 rpm		
	max. RCF	4427 xg	4427 xg	4427 xg		
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 12x15ml RB/FA, ID		15ml	17mm	RB/FA	12
1207353	Adapter, 1x1.5/2.0ml D11mm, 6/pk		1.5/2.0ml	11mm		12
1207345	Adapter, 1x7ml D13.5mm RB, 2/pk		7ml	13.5mm	RB	12
1207346	Adapter, 1x5ml D13.5mm, 2/pk		5ml	13.5mm		12


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo	Camlab product Code	1207259	Description	Rotor, Angle, 44x1.5/2.0ml, ID		
	Compatible with:		FC5718	FC5718R	FC5816	FC5816R
	Camlab product Code		1206876	1206877	1206878	1206879
	max. speed		13500 rpm	13500 rpm	15000 rpm	16000 rpm
	max. RCF		17113 xg	17113 xg	21128 xg	24039 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 44x1.5/2.0ml, ID		1.5/2.0ml	11mm		44
1207351	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		44
1207352	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		44


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

Rotor Photo	Camlab product Code	1207260	Description	Rotor, Angle, 30x1.5/2.0ml, ID, Sealable		
	Compatible with:	FC5714	FC5718	FC5718R		
	Camlab product Code	1206875	1206876	1206877		
	max. speed	12000 rpm	13000 rpm	14000 rpm		
	max. RCF	15131 xg	17758 xg	20595 xg		
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 30x1.5/2.0ml, ID, Sealable		1.5/2.0ml	11mm		30
1207351	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		30
1207352	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		30


Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

Rotor Photo 	Camlab product Code	1207261	Description	Rotor, Angle, 24x1.5/2.0ml, ID		
	Compatible with:	FC5714	FC5718	FC5718R	FC5816	FC5816R
	Camlab Product Code	1206875	1206876	1206877	1206878	1206879
	max. speed	14000 rpm	14000 rpm	15000 rpm	15000 rpm	16000 rpm
	max. RCF	18624 xg	18624 xg	21379 xg	21379 xg	24325 xg
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No.of tubes/rotor
No Adapter	Rotor, Angle, 24x1.5/2.0ml, ID		1.5/2.0ml	11mm		24
1207351	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		24
1207352	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		24


Rotors and Accessories for Ohaus Frontier FC5718, FC5718R

Rotor Photo 	Camlab product Code	1207262	Description	Rotor, Angle, 12x1.5/2.0ml, ID		
	Compatible with:		FC5718	FC5718R		
	Camlab product Code		1206876	1206877		
	max. speed		18000 rpm	18000 rpm		
	max. RCF		23542 xg	23542 xg		
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 12x1.5/2.0ml, ID		1.5/2.0ml	11mm		12
1207351	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		12
1207352	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		12
1207365	Adapter, 1x1.5ml for over16000xg, 6/pk		1.5ml	11mm		12

Rotors and Accessories for Ohaus Frontier FC5718, FC5718R

Rotor Photo 	Camlab product Code 1207263	Description	Rotor, Angle, 64x0.5ml, ID			
	Compatible with:		FC5718	FC5718R		
	Camlab product Code max. speed		1206876	1206877		
	max. RCF		13500 rpm	13500 rpm		
			16298 xg	16298 xg		
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 64x0.5ml, ID		0.5ml	8mm		64

Rotors and Accessories for Ohaus Frontier FC5718, FC5718R

	Camlab product Code	1207265	Description	Rotor, Angle, 4x8-w PCR Strip, ID		
	Compatible with:		FC5718	FC5718R		
	Camlab product Code		1206876	1206877		
	max. speed		15000 rpm	15000 rpm		
	max. RCF		15343 xg	15343 xg		
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 4x8-w PCR Strip, ID		8-w PCR Strip			32