




## Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5816, FC5816R

<b>Rotor Photo</b> 	<b>Camlab product Code</b>	<b>1207248</b>	Description	Rotor, Angle, 24x1.5/2.0ml, ID, Sealable		
	Compatible with:	FC5714	FC5718	FC5718R	FC5816	FC5816R
	<b>Camlab product Code</b>	<b>1206875</b>	<b>1206876</b>	<b>1206877</b>	<b>1206878</b>	<b>1206879</b>
	max. speed	14000 rpm	14000 rpm	15000 rpm	15000 rpm	16000 rpm
	max. RCF	18624 xg	18624 xg	21379 xg	21379 xg	24325 xg
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No.of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 24x1.5/2.0ml, ID, Sealable		1.5/2.0ml	11mm		24
<b>1207351</b>	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		24
<b>1207352</b>	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		24


## Rotors and Accessories for Ohaus Frontier FC5816, FC5816R

	<b>Camlab product Code</b>	<b>1207267</b>	Description	Rotor, Swing out, 4x250ml, ID		
	Compatible with:				FC5816	FC5816R
	<b>Camlab product Code</b>				<b>1206878</b>	<b>1206879</b>
	max. speed				4500 rpm	4500 rpm
	max. RCF				3780 xg	3780 xg
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>1207338</b>	Rack, 8x4-10ml D16mm Vac, 2/pk		4-10ml	16mm	Vac	32
<b>1207366</b>	Adapter, 9x1.5/2.0ml D11mm, 2/pk		1.5/2.0ml	11mm		36
<b>1207373</b>	Screw Cap, Rack of 250ml Rotor, 2/pk					
<b>1207374</b>	O-ring, Rack of 250ml Rotor, 10/pk					
<b>1207323</b>	Rack, 1x250ml D62mm FB, 2/pk		250ml	62mm	FB	4
<b>1207324</b>	Rack, 1x100ml D41mm RB, 2/pk		100ml	41mm	RB	4
<b>1207325</b>	Rack, 3x50ml D29mm RB, 2/pk		50ml	29mm	RB	12
<b>1207326</b>	Rack, 3x50ml D29mm FA, 2/pk		50ml	29mm	FA	12
<b>1207327</b>	Rack, 1x50ml D34mm RB, 2/pk, v3		50ml	34mm	RB	4
<b>1207328</b>	Rack, 2x50ml D29mm Rim, 2/pk		50ml	29mm	Rim	8
<b>1207329</b>	Rack, 1x50ml D29mm Rim, 2/pk, v3		50ml	29mm	Rim	4
<b>1207330</b>	Rack, 5x25ml D24.5mm RB, 2/pk		25ml	24.5mm	RB	20
<b>1207331</b>	Rack, 3x25ml D24.5mm Rim, 2/pk		25ml	24.5mm	Rim	12
<b>1207332</b>	Rack, 9x15ml D17mm RB, 2/pk		15ml	17mm	RB	36
<b>1207333</b>	Rack, 7x15ml D17mm FA, 2/pk		15ml	17mm	FA	28
<b>1207334</b>	Rack, 7x15ml D17mm Rim, 2/pk		15ml	17mm	Rim	28
<b>1207335</b>	Rack, 8x10ml D16mm RB, 2/pk		10ml	16mm	RB	32
<b>1207336</b>	Rack, 14x5-7ml D13mm RB, 2/pk		5-7ml	13mm	RB	56
<b>1207337</b>	Rack, 10x1.6-7ml D13mm Vac, 2/pk		1.6-7ml	13mm	Vac	40


## Rotors and Accessories for Ohaus Frontier FC5816, FC5816R

<b>Rotor Photo</b> 	<b>Camlab product Code</b>	<b>1207249</b>	Description	Rotor, Angle, 6x250ml FB, ID		
	Compatible with:				FC5816	FC5816R
	<b>Camlab product Code</b>				<b>1206878</b>	<b>1206879</b>
	max. speed				8000 rpm	8000 rpm
	max. RCF				10016 xg	10016 xg
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 6x250ml FB, ID		250ml	62mm	FB	6
<b>1207367</b>	Adapter, 8x1.5ml D11mm, 2/pk		1.5ml	11mm		48
<b>1207368</b>	Adapter, 5x10ml D16mm RB, 2/pk		10ml	16mm	RB	30
<b>1207369</b>	Adapter, 4x15ml D17mm FA, 2/pk		15ml	17mm	FA	24
<b>1207370</b>	Adapter, 2x30ml D26mm RB, 2/pk		30ml	26mm	RB	12
<b>1207371</b>	Adapter, 1x50ml D29mm RB, 2/pk		50ml	29mm	RB	6
<b>1207372</b>	Adapter, 1x50ml D30mm FA, 2/pk		50ml	30mm	FA	6


## Rotors and Accessories for Ohaus Frontier FC 5714, FC5816, FC5816R

	<b>Camlab product Code</b>	<b>1207268</b>	Description	Rotor, Swing out, 4x100ml, ID		
	Compatible with:	FC5714	FC5718	FC5718R		
	<b>Camlab product Code</b>	<b>1206875</b>	<b>1206876</b>	<b>1206877</b>		
	max. speed	4500 rpm	4500 rpm	4500 rpm		
	max. RCF	3350 xg	3350 xg	3350 xg		
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>1207289</b>	Rack, 1x100ml D46mm RB, 2/pk		100ml	46mm	RB	4
<b>1207290</b>	Rack, 1x85ml D38mm RB, 2/pk		85ml	38mm	RB	4
<b>1207291</b>	Rack, 1x50ml D34mm RB, 2/pk		50ml	34mm	RB	4
<b>1207292</b>	Rack, 1x50ml D29mm RB, 2/pk		50ml	29mm	RB	4
<b>1207293</b>	Rack, 1x50ml D29mm FA, 2/pk		50ml	29mm	FA	4
<b>1207294</b>	Rack, 1x50ml D29mm Rim, 2/pk		50ml	29mm	Rim	4
<b>1207295</b>	Rack, 1x30ml D25mm Rim, 2/pk		30ml	25mm	Rim	4
<b>1207296</b>	Rack, 1x25ml D25mm RB, 2/pk		25ml	25mm	RB	4
<b>1207297</b>	Rack, 2x15ml D17mm FA, 2/pk		15ml	17mm	FA	8
<b>1207298</b>	Rack, 1x15ml D17mm Rim, 2/pk		15ml	17mm	Rim	4
<b>1207299</b>	Rack, 7x15ml D17mm RB, 2/pk		15ml	17mm	RB	28
<b>1207300</b>	Rack, 9x5-7ml D13mm RB, 2/pk		5-7ml	13mm	RB	36
<b>1207301</b>	Rack, 10x1.5ml D11mm, 2/pk		1.5ml	11mm		40
<b>1207302</b>	Rack, 9x4.5ml D12mm FB, 2/pk		4.5ml	12mm	FB	36
<b>1207303</b>	Rack, 7x5-7ml D13mm RB, 2/pk		5-7ml	13mm	RB	28
<b>1207304</b>	Rack, 7x5-7ml D13mm Vac, 2/pk		5-7ml	13mm	Vac	28
<b>1207305</b>	Rack, 5x4-10ml D16mm Sar, 2/pk		4-10ml	16mm	Sar	20


## Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

<b>Rotor Photo</b> 	<b>Camlab product Code</b>	<b>1207269</b>	Description	Rotor, Swing out, 4x100ml, ID, Sealable		
	Compatible with:	<b>FC5714</b>	<b>FC5718</b>	<b>FC5718R</b>		
	<b>Camlab product Code</b>	<b>1206875</b>	<b>1206876</b>	<b>1207877</b>		
	max. speed	4000 rpm	5000 rpm	5000 rpm		
	max. RCF	2486 xg	3885 xg	3885 xg		
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>1207283</b>	Bucket, 100ml w/o Cap, 2/pk					
<b>1207284</b>	Bucket, 100ml w/ Cap, Sealable, 2/pk					
<b>1207375</b>	Screw Cap, 100ml Bucket, Sealable, 2/pk					
<b>1207376</b>	O-ring, 100ml Bucket, Sealable, 10/pk					
<b>1207306</b>	Rack, 1x100ml D40mm RB, 2/pk		100ml	40mm	RB	4
<b>1207307</b>	Rack, 1x85ml D38mm RB, 2/pk, v2		85ml	38mm	RB	4
<b>1207308</b>	Rack, 1x50ml D34mm RB, 2/pk, v2		50ml	34mm	RB	4
<b>1207309</b>	Rack, 1x50ml D29mm RB, 2/pk, v2		50ml	29mm	RB	4
<b>1207310</b>	Rack, 1x50ml D29mm FA, 2/pk, v2		50ml	29mm	FA	4
<b>1207311</b>	Rack, 1x50ml D29mm Rim, 2/pk, v2		50ml	29mm	Rim	4
<b>1207312</b>	Rack, 1x30ml D25mm Rim, 2/pk, v2		30ml	25mm	Rim	4
<b>1207313</b>	Rack, 1x25ml D25mm RB, 2/pk, v2		25ml	25mm	RB	4
<b>1207314</b>	Rack, 1x15ml D17mm FA, 2/pk		15ml	17mm	FA	4
<b>1207315</b>	Rack, 2x15ml D17mm FA, 2/pk, v2		15ml	17mm	FA	8
<b>1207316</b>	Rack, 2x15ml D17mm Rim, 2/pk		15ml	17mm	Rim	8
<b>1207317</b>	Rack, 4x15ml D17mm RB, 2/pk		15ml	17mm	RB	16
<b>1207318</b>	Rack, 7x5-7ml D13mm RB, 2/pk, v2		5-7ml	13mm	RB	28
<b>1207319</b>	Rack, 5x1.5ml D11mm, 2/pk		1.5ml	11mm		20
<b>1207320</b>	Rack, 3x10ml D17mm Sar, 2/pk		10ml	17mm	Sar	12
<b>1207321</b>	Rack, 4x5-7ml D13mm RB, 2/pk		5-7ml	13mm	RB	16
<b>1207322</b>	Rack, 4x10ml D17mm Vac, 2/pk		10ml	17mm	Vac	16
<b>1207285</b>	Bucket, 2x50ml D29mm FA, w/ Rack, 2/pk		50ml	29mm	FA	8
<b>1207286</b>	Bucket, 3x15ml D17mm FA, w/ Rack, 2/pk		15ml	17mm	FA	12
<b>1207287</b>	Cyto, Rectangular Bucket, 2/pk					
<b>1207377</b>	Cyto, Slide Carrier, 1/pk					
<b>1207378</b>	Cyto, Container, 2x1.5ml Cyl., 1/pk					
<b>1207379</b>	Cyto, Container, 2x3ml Conical, 1/pk					
<b>1207380</b>	Cyto, Container Cover, 1/pk					
<b>1207381</b>	Cyto, Filter Card, Cyl., 100/pk					
<b>1207382</b>	Cyto, Filter Card, Conical, 100/pk					


## Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

<b>Rotor Photo</b> 	<b>Camlab product Code</b>	<b>1207270</b>	Description	Rotor, Swing out, 2x3MTP, ID		
	Compatible with:	<b>FC5714</b>	<b>FC5718</b>	<b>FC5718R</b>	<b>FC5816</b>	<b>FC5816R</b>
	<b>Camlab Product Code</b>	<b>1206875</b>	<b>1206876</b>	<b>1207877</b>	<b>1206878</b>	<b>1206879</b>
	max. speed	4500 rpm	4500 rpm	4500 rpm	4500 rpm	4500 rpm
	max. RCF	2716 xg	2716 xg	2716 xg	2716 xg	2716 xg
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Swing out, 2x3MTP, ID		MTP			6
<b>1207288</b>	Bucket, 3xMTP, 2/pk		MTP			6
<b>1207356</b>	Adapter, 48x1.5/2.0ml D11mm, 2/pk		1.5/2.0ml	11mm		96

## Rotors and Accessories for Ohaus Frontier FC5718, FC5718


Rotor Photo	Camlab product Code	1207250	Description	Rotor, Angle, 6x85ml RB, ID, Hi		
	Compatible with		<b>FC5718</b>	<b>FC5718R</b>		
	Camlab product Code		<b>1206876</b>	<b>1206877</b>		
	max. speed		11000 rpm	13500 rpm		
	max. RCF		13932 xg	20984 xg		
Accessories						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
No Adapter	Rotor, Angle, 6x85ml RB, ID, Hi		85ml	38mm	RB	6
1207357	Adapter, 1x15ml D17mm RB, 2/pk		15ml	17mm	RB	6
1207358	Adapter, 1x30ml D25mm RB, 2/pk		30ml	25mm	RB	6
1207359	Adapter, 1x50ml D29mm RB, 2/pk, v2		50ml	29mm	RB	6
1207360	Adapter, 1x50ml D29mm FA, 2/pk		50ml	29mm	FA	6
1207361	Adapter, 1x15ml D17mm FA, 2/pk		15ml	17mm	FA	6
1207362	Adapter, 2x15ml D17mm RB, 2/pk		15ml	17mm	RB	12
1207363	Adapter, 1x16ml D18mm RB, 2/pk		16ml	18mm	RB	6

## Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R


<b>Rotor Photo</b> 	<b>Camlab product Code</b> <b>1207251</b>	Description	Rotor, Angle, 6x85ml RB, ID			
	Compatible with:	<b>FC5718</b>	<b>FC5718R</b>	<b>FC5816</b>	<b>FC5816R</b>	
	<b>Camlab product Code</b> max. speed	<b>1206876</b>	<b>1206877</b>	<b>1206878</b>	<b>1206879</b>	
	max. RCF	9000 rpm	9000 rpm	11000 rpm	13000 rpm	
		10413 xg	10413 xg	15555 xg	21726 xg	
<b>Accessories</b>						
<b>Camlab product Code</b>	Description	Vol./Tube	Diameter	Tube type	No.of tubes/rotor	
<b>No Adapter</b>	Rotor, Angle, 6x85ml RB, ID	85ml	38mm	RB	6	
<b>1207357</b>	Adapter, 1x15ml D17mm RB, 2/pk	15ml	17mm	RB	6	
<b>1207358</b>	Adapter, 1x30ml D25mm RB, 2/pk	30ml	25mm	RB	6	
<b>1207359</b>	Adapter, 1x50ml D29mm RB, 2/pk, v2	50ml	29mm	RB	6	
<b>1207360</b>	Adapter, 1x50ml D29mm FA, 2/pk	50ml	29mm	FA	6	
<b>1207361</b>	Adapter, 1x15ml D17mm FA, 2/pk	15ml	17mm	FA	6	
<b>1207362</b>	Adapter, 2x15ml D17mm RB, 2/pk	15ml	17mm	RB	12	
<b>1207363</b>	Adapter, 1x16ml D18mm RB, 2/pk	16ml	18mm	RB	6	




## Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo	Camlab product Code	1207252	Description	Rotor, Angle, 4x85ml RB, ID, Hi		
	Compatible with:		<b>FC5718</b>	<b>FC5718R</b>	<b>FC5816</b>	<b>FC5816R</b>
	Camlab product Code		<b>1206876</b>	<b>1206877</b>	<b>1206878</b>	<b>1206879</b>
	max. speed		12000 rpm	12000 rpm	12000 rpm	12000 rpm
	max. RCF		14809 xg	14809 xg	14809 xg	14809 xg
<b>Accessories</b>						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 4x85ml RB, ID, Hi		85ml	38mm	RB	4
<b>1207357</b>	Adapter, 1x15ml D17mm RB, 2/pk		15ml	17mm	RB	4
<b>1207358</b>	Adapter, 1x30ml D25mm RB, 2/pk		30ml	25mm	RB	4
<b>1207359</b>	Adapter, 1x50ml D29mm RB, 2/pk, v2		50ml	29mm	RB	4
<b>1207360</b>	Adapter, 1x50ml D29mm FA, 2/pk		50ml	29mm	FA	4
<b>1207361</b>	Adapter, 1x15ml D17mm FA, 2/pk		15ml	17mm	FA	4
<b>1207362</b>	Adapter, 2x15ml D17mm RB, 2/pk		15ml	17mm	RB	8
<b>1207363</b>	Adapter, 1x16ml D18mm RB, 2/pk		16ml	18mm	RB	4


## Rotors and Accessories for Ohaus Frontier FC5816, FC5816R

<b>Rotor Photo</b> 	<b>Camlab product Code</b>	<b>1207271</b>	Description	Rotor, Swing out, 16x50ml, ID		
	Compatible with:				<b>FC5816</b>	<b>FC5816R</b>
	<b>Camlab product Code</b>				<b>1206878</b>	<b>1206879</b>
	max. speed				4500 rpm	4500 rpm
	max. RCF				3735 xg	3735 xg
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>1207339</b>	Rack, 4x50ml D29mm FA, 2/pk		50ml	29mm	FA	16
<b>1207340</b>	Rack, 10x15ml D17mm FA, 2/pk		15ml	17mm	FA	40
<b>1207341</b>	Rack, 14x5-7ml D13mm Vac/Sar, 2/pk		5-7ml	13mm	Vac/Sar	56
<b>1207342</b>	Rack, 12x4-10ml D16mm Vac/Sar, 2/pk		4-10ml	16mm	Vac/Sar	48


## Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo	Camlab product Code	1207253	Description	Rotor, Angle, 10x50ml FA, ID		
	Compatible with:		<b>FC5718</b>	<b>FC5718R</b>	<b>FC5816</b>	<b>FC5816R</b>
	Camlab product Code		<b>1206876</b>	<b>1206877</b>	<b>1206878</b>	<b>1206879</b>
	max. speed		7500 rpm	7500 rpm	9000 rpm	10500 rpm
	max. RCF		8174 xg	8174 xg	11771 xg	16022 xg
<b>Accessories</b>						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 10x50ml FA, ID		50ml	29mm	FA	10
<b>1207353</b>	Adapter, 1x1.5/2.0ml D11mm, 6/pk		1.5/2.0ml	11mm		10
<b>1207345</b>	Adapter, 1x7ml D13.5mm RB, 2/pk		7ml	13.5mm	RB	10
<b>1207346</b>	Adapter, 1x5ml D13.5mm, 2/pk		5ml	13.5mm		10
<b>1207350</b>	Adapter, 1x15ml D17mm RB, 2/pk, v2		15ml	17mm	RB	10


## Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

<b>Rotor Photo</b> 	<b>Camlab product Code</b>	<b>1207254</b>	Description	Rotor, Angle, 6x50ml RB/FA, ID		
	Compatible with:	<b>FC5714</b>	<b>FC5718</b>	<b>FC5718R</b>		
	<b>Camlab product Code</b>	<b>1206875</b>	<b>1206876</b>	<b>1206877</b>		
	max. speed	6000 rpm	6000 rpm	6000 rpm		
	max. RCF	4427 xg	4427 xg	4427 xg		
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 6x50ml RB/FA, ID		50ml	29mm	RB/FA	6
<b>1207353</b>	Adapter, 1x1.5/2.0ml D11mm, 6/pk		1.5/2.0ml	11mm		6
<b>1207345</b>	Adapter, 1x7ml D13.5mm RB, 2/pk		7ml	13.5mm	RB	6
<b>1207346</b>	Adapter, 1x5ml D13.5mm, 2/pk		5ml	13.5mm		6
<b>1207347</b>	Adapter, 1x30ml D26mm, 2/pk		30ml	26mm		6
<b>1207348</b>	Adapter, 1x16ml D18mm, 2/pk		16ml	18mm		6
<b>1207349</b>	Adapter, 1x15ml D17mm, 2/pk		15ml	17mm		6
<b>1207350</b>	Adapter, 1x15ml D17mm RB, 2/pk, v2		15ml	17mm	RB	6


## Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

<b>Rotor Photo</b> 	<b>Camlab product Code</b>	<b>1207255</b>	Description	Rotor, Angle, 6x50ml RB, ID, Hi		
	Compatible with		<b>FC5718</b>	<b>FC5718R</b>	<b>FC5816</b>	<b>FC5816R</b>
	<b>Camlab product Code</b>		<b>1206876</b>	<b>1206877</b>	<b>1206878</b>	<b>1206879</b>
	max. speed		12000 rpm	12000 rpm	13000 rpm	13000 rpm
	max. RCF		13522 xg	13522 xg	15869 xg	15869 xg
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 6x50ml RB, ID, Hi		50ml	29mm	RB	6
<b>1207347</b>	Adapter, 1x30ml D26mm, 2/pk		30ml	26mm		6
<b>1207348</b>	Adapter, 1x16ml D18mm, 2/pk		16ml	18mm		6
<b>1207349</b>	Adapter, 1x15ml D17mm, 2/pk		15ml	17mm		6
<b>1207364</b>	Adapter, 1x15ml D17.5mm FA, 2/pk		15ml	17.5mm	FA	6


## Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

<b>Rotor Photo</b> 	<b>Camlab product Code</b>	<b>1207256</b>	Description	Rotor, Angle, 30x15ml RB/FA, ID		
	Compatible with:	<b>FC5714</b>	<b>FC5718</b>	<b>FC5718R</b>	<b>FC5816</b>	<b>FC5816R</b>
	<b>Camlab Product Code</b>	<b>1206875</b>	<b>1206876</b>	<b>1206877</b>	<b>1206878</b>	<b>1206879</b>
	max. speed	4500 rpm	4500 rpm	4500 rpm	4500 rpm	4500 rpm
	max. RCF	2830 xg	2830 xg	2830 xg	2830 xg	2830 xg
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No.of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 30x15ml RB/FA, ID		15ml	17mm	RB/FA	30
<b>1207353</b>	Adapter, 1x1.5/2.0ml D11mm, 6/pk		1.5/2.0ml	11mm		30
<b>1207345</b>	Adapter, 1x7ml D13.5mm RB, 2/pk		7ml	13.5mm	RB	30
<b>1207346</b>	Adapter, 1x5ml D13.5mm, 2/pk		5ml	13.5mm		30

## Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R


<b>Rotor Photo</b> 	<b>Camlab product Code</b>	<b>1207257</b>	Description	Rotor, Angle, 20x10ml RB, ID, Hi		
	Compatible with:		<b>FC5718</b>	<b>FC5718R</b>	<b>FC5816</b>	<b>FC5816R</b>
	<b>Camlab Product Code</b>		<b>1206876</b>	<b>1206877</b>	<b>1206878</b>	<b>1206879</b>
	max. speed		12000 rpm	12000 rpm	12000 rpm	12000 rpm
	max. RCF		15775 xg	15775 xg	15775 xg	15775 xg
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 20x10ml RB, ID, Hi		10ml	16mm	RB	20

## Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R


<b>Rotor Photo</b> 	<b>Camlab product Code</b>	<b>1207258</b>	Description	Rotor, Angle, 12x15ml RB/FA, ID		
	Compatible with:	<b>FC5714</b>	<b>FC5718</b>	<b>FC5718R</b>		
	<b>Camlab Product Code</b>	<b>1206875</b>	<b>1206876</b>	<b>1206877</b>		
	max. speed	6000 rpm	6000 rpm	6000 rpm		
	max. RCF	4427 xg	4427 xg	4427 xg		
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 12x15ml RB/FA, ID		15ml	17mm	RB/FA	12
<b>1207353</b>	Adapter, 1x1.5/2.0ml D11mm, 6/pk		1.5/2.0ml	11mm		12
<b>1207345</b>	Adapter, 1x7ml D13.5mm RB, 2/pk		7ml	13.5mm	RB	12
<b>1207346</b>	Adapter, 1x5ml D13.5mm, 2/pk		5ml	13.5mm		12




## Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816, FC5816R

Rotor Photo	Camlab product Code	1207259	Description	Rotor, Angle, 44x1.5/2.0ml, ID		
	Compatible with:		<b>FC5718</b>	<b>FC5718R</b>	<b>FC5816</b>	<b>FC5816R</b>
	Camlab product Code		<b>1206876</b>	<b>1206877</b>	<b>1206878</b>	<b>1206879</b>
	max. speed		13500 rpm	13500 rpm	15000 rpm	16000 rpm
	max. RCF		17113 xg	17113 xg	21128 xg	24039 xg
<b>Accessories</b>						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 44x1.5/2.0ml, ID		1.5/2.0ml	11mm		44
<b>1207351</b>	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		44
<b>1207352</b>	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		44


## Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R

Rotor Photo	Camlab product Code	1207260	Description	Rotor, Angle, 30x1.5/2.0ml, ID, Sealable		
	Compatible with:	<b>FC5714</b>	<b>FC5718</b>	<b>FC5718R</b>		
	Camlab product Code	<b>1206875</b>	<b>1206876</b>	<b>1206877</b>		
	max. speed	12000 rpm	13000 rpm	14000 rpm		
	max. RCF	15131 xg	17758 xg	20595 xg		
<b>Accessories</b>						
Camlab product Code	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 30x1.5/2.0ml, ID, Sealable		1.5/2.0ml	11mm		30
<b>1207351</b>	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		30
<b>1207352</b>	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		30


## Rotors and Accessories for Ohaus Frontier FC5714, FC5718, FC5718R, FC5816, FC5816R

<b>Rotor Photo</b> 	<b>Camlab product Code</b>	1207261	Description	Rotor, Angle, 24x1.5/2.0ml, ID		
	Compatible with:	<b>FC5714</b>	<b>FC5718</b>	<b>FC5718R</b>	<b>FC5816</b>	<b>FC5816R</b>
	<b>Camlab Product Code</b>	<b>1206875</b>	<b>1206876</b>	<b>1206877</b>	<b>1206878</b>	<b>1206879</b>
	max. speed	14000 rpm	14000 rpm	15000 rpm	15000 rpm	16000 rpm
	max. RCF	18624 xg	18624 xg	21379 xg	21379 xg	24325 xg
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No.of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 24x1.5/2.0ml, ID		1.5/2.0ml	11mm		24
<b>1207351</b>	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		24
<b>1207352</b>	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		24


## Rotors and Accessories for Ohaus Frontier FC5718, FC5718R

<b>Rotor Photo</b> 	<b>Camlab product Code</b>	1207262	Description	Rotor, Angle, 12x1.5/2.0ml, ID		
	Compatible with:		<b>FC5718</b>	<b>FC5718R</b>		
	<b>Camlab product Code</b>		<b>1206876</b>	<b>1206877</b>		
	max. speed		18000 rpm	18000 rpm		
	max. RCF		23542 xg	23542 xg		
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 12x1.5/2.0ml, ID		1.5/2.0ml	11mm		12
<b>1207351</b>	Adapter, 1x0.2/0.4ml D6mm, 6/pk		0.2/0.4ml	6mm		12
<b>1207352</b>	Adapter, 1x0.5ml D8mm, 6/pk		0.5ml	8mm		12
<b>1207365</b>	Adapter, 1x1.5ml for over16000xg, 6/pk		1.5ml	11mm		12

### Rotors and Accessories for Ohaus Frontier FC5718, FC5718R, FC5816R

<b>Rotor Photo</b> 	<b>Camlab product Code</b> 1207263	Description	Rotor, Angle, 64x0.5ml, ID			
	Compatible with:		<b>FC5718</b>	<b>FC5718R</b>		
	<b>Camlab product Code</b> max. speed		<b>1206876</b> 13500 rpm	<b>1206877</b> 13500 rpm		
	max. RCF		16298 xg	16298 xg		
<b>Accessories</b>						
<b>Camlab product Code</b>	Description	Vol./Tube	Diameter	Tube type	No. of tubes/rotor	
<b>No Adapter</b>	Rotor, Angle, 64x0.5ml, ID	0.5ml	8mm		64	

## Rotors and Accessories for Ohaus Frontier FC5718, FC5718R

	<b>Camlab product Code</b>	<b>1207265</b>	Description	Rotor, Angle, 4x8-w PCR Strip, ID		
	Compatible with:		<b>FC5718</b>	<b>FC5718R</b>		
	<b>Camlab product Code</b>		<b>1206876</b>	<b>1206877</b>		
	max. speed		15000 rpm	15000 rpm		
	max. RCF		15343 xg	15343 xg		
<b>Accessories</b>						
<b>Camlab product Code</b>	Description		Vol./Tube	Diameter	Tube type	No. of tubes/rotor
<b>No Adapter</b>	Rotor, Angle, 4x8-w PCR Strip, ID		8-w PCR Strip			32