

# Osmometers at a glance



Enter this number into the search box at [www.camlab.co.uk](http://www.camlab.co.uk) to see more information, prices and live stock data

Our Barcode	Product	Sample Volume	Measuring Time	Reproducibility	Measuring Range Mosm / kg H <sub>2</sub> O	Freezing Operation
1138947	Micro-Digital Osmometer Type 6	50 – 100 µL	1.5 minutes	+/- 0.5% or 1.5 mosm	0 - 2500	Manual freezing initiation
1192570	Micro-Digital Osmometer Type 6M	10 – 50 µL	1.3 minutes	+/- 0.5% or 1.5 mosm with 50 uL sample  +/-1% or 3 mosm with 10-25 uL sample	0 - 2500	Manual freezing initiation
1138940	Micro-Digital Osmometer MOD200 Plus	50 – 100 µL	1.5 minutes	+/-0.5% or 1.5 mosm	0 - 2500	Automatic Process
1192573	Micro-Digital Osmometer Type 15M	10 – 50 µL	1.5 minutes	+/- 0.5% or 1.5 mosm with 50 uL sample  +/-1% or 3 mosm with 10-25 uL sample	0 - 2500	Automatic Process
V1.0	\\WEBSERVER\storesphere_folders\siteimages\0.Quick Guides Uploaded to Web					