

Brannan Product Data Sheet

PDS0501.20.REV1

SPH011 Hydrometer - Special Purpose Hydrometer - Latexometer



These instruments are calibrated in grams of dry rubber per litre of latex and adjusted at 28.9°C (84°F). 250mm long.

All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.

Product Features:

- Calibrated in grams of dry rubber per litre of latex
- Adjusted at 28.9°C (84°F)
- 250mm long
- Made in UK
- Supplied in protective, resealable plastic tube

Product Specifications:

Material:	Soda glass, steel shot ballast, synthetic wax resin
Accuracy:	+/- 10.0
Ranges:	50 to 450 g/l or 50 to 250 g/l
Bulb Diameter:	26mm or 43mm
Overall length:	250mm

Standard Product Numbers

Range	Bulb Diameter	Gross Product Weight*	Product Number
50 to 450 g/l	26mm	100g	10/611/0
50 to 250 g/l	43mm	133g	10/612/0

*including resealable plastic tube

www.brannan.co.uk

Measurement for Life.

