

# Brannan Product Data Sheet

PDS0517.20.REV1

## SPH012 Hydrometer - Special Purpose Hydrometer - Urinometer



These hydrometers are 160mm long and measure the specific gravity of urine at 15.6°C (60°F). All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.

### Product Features:

- Measures the specific gravity of urine
- Adjusted at 15.6°C (60°F)
- 160mm long
- Made in UK
- Supplied in protective, resealable plastic tube

### Product Specifications:

<b>Material:</b>	Soda glass, steel shot ballast, synthetic wax resin
<b>Range:</b>	1.000 to 1.060°SG
<b>Accuracy:</b>	+/- 0.002°SG
<b>Divisions:</b>	0.001°SG
<b>Bulb Diameter:</b>	16mm
<b>Overall Length:</b>	160mm
<b>Weight:</b>	21g (including protective plastic tube)

### Standard Product Numbers

Range	Divisions	Product No
1.000 to 1.060°SG	0.001°SG	10/650/0

[www.brannan.co.uk](http://www.brannan.co.uk)

**Measurement for Life.**

