



Funnels, dropping, graduated, PTFE key

MBL®

- Complies to ISO 4800
- Manufactured from borosilicate glass
- With double graduated scale showing approx. volume contained/volume dispensed
- Graduations in white enamel
- With interchangeable PTFE stopcock for grease-free performance
- Supplied with polyethylene stopper

Product Code	Capacity (ml)	Grads. ml	Lowest grad. mark ml	Stopcock bore mm	Stopper size	Pack qty.
FST511	50	1	5	1.5	19/26	1
FST521	100	2	10	1.5	19/26	1
FST531	250	5	25	2.5	24/29	1
FST542	500	10	50	4.0	24/29	1
FST552	1000	20	100	4.0	29/32	1



Funnels, dropping, graduated, PTFE key, with cone

PYREX®

- Complies with ISO 4800
- With interchangeable PTFE key for grease free performance
- Double graduated scale showing volume of liquid contained/delivered
- Stem inside cone for easy fitting to Quickfit® assemblies
- Graduations in white enamel

Product Code	Capacity (ml)	Grads. ml	Stopcock bore mm	Cone and socket size	Pack qty.
443333	50	1	1.5	19/26	1
443340	100	2	1.5	19/26	1
443347	250	5	2.5	29/32	1
443354	500	10	4.0	29/32	1
443360	1000	20	4.0	29/32	1