




Filtration system, for 47mm diameter membranes



- All glass filtration system for use with 47mm diameter membranes
- Borosilicate glass sintered disc to support membrane
- Comprises four components which should be ordered separately
 - Funnel 300ml capacity
 - Support base with sintered disc and vacuum connection
 - Metal retaining clip
 - 1000ml conical flask suitable for use under vacuum
- Fully autoclavable
- Ideal for microbiology



Item Code	Description	Pack qty.
5809/2	Funnel	1
5810/3	Base	1
5810/4	Flask	1
5810/5	Clip	1

Funnels, lightweight, polypropylene *Azlon*

- Lightweight polypropylene funnels with ribs that prevent air locks, allowing rapid filtration
- Autoclavable to 121°C
- Made from food grade plastic



Item Code	Top I.D. mm	Stem O.D. at tip mm	Overall height mm	Pack qty.
FWC052	30	2	47	10
FWC054	30	5	47	10
FWC056	40	5	62	10
FWC058	50	9	78	10
FWC060	75	10	107	10
FWC062	100	13	149	10
FWC064	120	14	170	5
FWC066	150	17	210	5

i General purpose plastic funnels available on page 175.
For information on food safety, please refer to the technical section.