



Tubes, air/steam inlet, 90° bend



- With detachable plastic screwthread connector to facilitate the attachment of 9mm bore flexible tubing
- Ground glass cone on stem fits Quickfit® sockets
- Lower end of tube left plain for Glassblower to cut to length and draw to jet

| Product Code | Cone size | Approx. stem length mm | Pack qty. |
|--------------|-----------|------------------------|-----------|
| MF15/1B/SC | 14/23 | 150 | 1 |
| MF15/2B/SC | 19/26 | 150 | 1 |
| MF15/3B/SC | 24/29 | 150 | 1 |



Tube, centrifuge, conical base, graduated

PYREX®

- Centrifuge tube with white enamel graduations
- Graduated 0 - 3ml in 0.1ml divisions and 3 - 10ml in 0.5ml divisions
- Manufactured from Pyrex® borosilicate glass
- Can be centrifuged up to 3000 RCF values

| Product Code | Capacity (ml) | Highest mark ml | Diam. mm | Length (mm) | Pack qty. |
|--------------|---------------|-----------------|----------|-------------|-----------|
| 3460/02M | 15 | 10 | 17 | 110 | 40 |



Always examine centrifuge tubes before use. Do not use if scratched, chipped or abraded.

Tube, centrifuge, conical base, plain

PYREX®

- Manufactured from chemically inert Pyrex® glass for contamination-free performance
- Conical shape with rim
- Can be centrifuged up to 3000 RCF values

| Product Code | Capacity (ml) | Diam. mm | Length (mm) | Pack qty. |
|--------------|---------------|----------|-------------|-----------|
| 1680/02M | 15 | 17 | 110 | 40 |

