




## Bottles, media-lab, without cap or pouring ring

**PYREX®**

- Glass bottles only
- As Pyrex 1516: range but supplied without cap or pouring ring
- Printed with retrace code: downloadable batch certificate

Item Code	Capacity ml	Bottle height mm	Body diam. mm	Thread size mm	Pack qty.
1517/01	25	70	36	25	10
1517/02	50	88	46	32	10
1517/04	100	100	56	45	10
1517/06	250	138	70	45	10
1517/08	500	176	86	45	10
1517/10	1000	225	101	45	10
1517/12	2000	260	136	45	10
1517/16	5000	335	181	45	1
1517/18	10000	410	227	45	1
1517/22	20000	510	300	45	1

 Choice of caps and pouring rings available separately on pages 56-57.

 Never autoclave bottles with a fully tightened screwcap. Always slacken off. Failure to follow this procedure can lead to breakage of the bottle.




## Bottles, media-lab, plastic-coated

**PYREX®**

- Bottles only - supplied without cap or pouring ring
- Borosilicate glass with clear plastic coating to provide added protection against mechanical breakage
- Reduces the hazard of flying glass and liquid splashes in the event of breakage
- Printed with retrace code: downloadable batch certificate
- Coating withstands steam sterilisation up to 135°C

Item Code	Capacity ml	Bottle height mm	Body diam. mm	Thread size mm	Pack qty.
1518/04	100	100	56	45	10
1518/06	250	138	70	45	10
1518/08	500	176	86	45	10
1518/10	1000	225	101	45	10
1518/12	2000	260	136	45	10
1518/16	5000	335	181	45	1
1518/18	10000	410	227	45	1

 For further information about plastic coated bottles please see the technical information section.

