



Evaluation Software

FOR PROGRAMMING AND READOUT OF EBRO DATA LOGGERS

-ebro-
a xylem brand

Evaluation Software

for ebro Data Loggers

For programming and readout of ebro data loggers and for evaluating the measurement values ebro offers three different software versions: the **free Winlog.basic**, the **standard software Winlog.light** and the **professional software Winlog.pro**.

Benefits


- Easy installation
- Easy programming of the data loggers, no prior knowledge required
- Extensive and custom report generation
- Suitable for all applications
- Security by compliance with FDA 21 CFR Part 11

Extensive Report Generation

The software makes it easy to generate standard and custom reports:


- Compact, one-sided report (1)
- Multi-page, detailed report (2)
- Tabular report with the measurement values (3)
- Insertion of your own company logo possible (4)
- Export data to Microsoft Excel® and PDF (5)
- Integration of pictures and graphics possible (Winlog.pro only) (6)

(1)




(2)

(3)



(4)

(5)



(6)



Of the Software Supported Data Logger Families

The table below shows the software version which can be used to evaluate the individual data logger family:

Supported data logger families	Winlog.basic	Winlog.light	Winlog.pro
EBI 10		✓	✓
EBI 11		✓	✓
EBI 15			✓
EBI 20	✓	✓	✓
EBI 30			✓
EBI 40		✓	✓
EBI 100		✓	✓
EBI 300 / EBI 310 / EBI 330	✓	✓	✓

Localized Languages

The software version Winlog.basic, Winlog.light and Winlog.pro are available in the following languages:

- German
 - English
 - French
- Italian
 - Spanish
 - Chinese
- Japanese
 - Czech
 - Swedish

The Different Software Versions at a Glance



Winlog.basic

The **easy to use, free software**

- User friendly: intuitive and easy to use
- Graphical and numerical display of measurement data
- Protocol print (with printer and computers)
- Scan and enlarge of the measurement data
- Data export to Microsoft® Excel and PDF
- 21 CFR Part 11 - functionality as an option



Winlog.light

The **standard software**

- Contains all the features of Winlog.basic
- Data import from other ebro application programs
- Many standard reports
- Displays statistics on measurement data (e.g. MIN/MAX, mean, standard deviation etc.)
- Creation of configurations possible

How to order:

Type	Description	Part No.
Winlog.basic	Free evaluation software	1340-2375
Winlog.light	Standard evaluation software	1340-2354
Winlog.pro	Professional evaluation software	1340-2355


+49 841 95478-0



Winlog.pro

The **professional software**

- Contains all the features of Winlog.basic and Winlog.light
- Enables real-time monitoring with wireless data loggers
- Formula editor for calculating the F0-value of the absolute humidity, the PE value etc.
- Display of the timeline either absolute or relative
- Customized definition of individual areas possible (with their own statistics and calculations)
- Including calibration tool for data loggers
- Integration of pictures and graphics into reports possible
- IQ/OQ documentation optional



System Requirements

So that the software can run on your computer with any problem, your computer must meet the following requirements:

Hardware Requirements:

- At least 1 GHz processor speed
- At least 1 GB working memory
- At least 300 MB free hard disc space
- USB (Universal Serial Bus)

Software requirements:

Operating system Microsoft®

- Windows XP
- Windows Vista (32 Bit and 64 Bit)
- Windows 7 (32 Bit and 64 Bit)



Distributed and Supported in the UK by:



24 Norman Way Industrial Estate, Over, Cambridge, CB24 5WE
Tel: +44(0)1954 233 100 Fax: +44(0)1954 233 101
Email: sales@camlab.co.uk Web: www.camlab.co.uk

Distributed by