

## Data Sheet

### UVSL-26P - 6-Watt, 254/365nm Rechargeable UV Lamp

#### Part Numbers:

95-0181-01 (115V)	- Non Ce Certified
95-0181-02 (230V)	- CE Certified

#### Intensity:

**Shortwave: 1 Tube** - Approximately 1130  $\mu\text{W}/\text{cm}^2$  at a distance of 3" from the sample

**Shortwave: 2 Tubes** - Approximately 1250  $\mu\text{W}/\text{cm}^2$  at a distance of 3" from the sample

#### Dimensions:

140mm Long x 76mm Wide x 249mm High

#### Weight:

1.8Kg

#### Spare Parts:

34-0034-01	- 6-Watt, 365nm UV Tube
34-0013-01	- 6-Watt, 254nm UV Tube

#### Accessories:

58-0127-02	- 12V Adapter
------------	---------------