



VP 1HV

Portable Bench Top Vacuum Pump



Twin Head Diaphragm Vacuum Pump

- For clean evacuation and transfer of air and gases.
- High performance and long working life through the stress-optimised KNF diaphragm.

Maximum Permissible gas and ambient temp: 40°C
Standard motor: 230V/50Hz

Features

- Contamination free pumping of air and gases
- Quiet and smooth running
- High tolerance of both vapour and condensate
- Simple operation
- Completely portable
- Maintenance-free

Model Codes and Materials

Type	Pump Head	Diaphragm	Valves
VP1HV	PPS	EPDM	EPDM

Performance Data

Type	Flow Rate Litres/min	Vacuum Mbar/abs	Pressure Barg
VP1HV	5	25	0.3

Motor Data

Type	Voltage/ Frequency	Power PI (W)	Operating Current (A)
VP1HV	230/50	65	0.65

Dimensions

Type	Weight	Dimensions
VP1HV	2.0 Kg	185 x 125 x 145 mm

Accessories

Description	Order Number
Service Kit for VP1HV	032498

