

BOD/Respiration

Distributed By:



**Dilution BOD,
BOD self-check measurement,
Respiration, Biogas determination**

| Application Range | | |
|-----------------------|---------------|-----------------------------------------------------------------------------------------------------------------------|
| | OxiTop®-i | OxiTop®-IDS with Multi 3620 IDS/Multi 3630 IDS |
| Application | BOD routine | BOD routine, standard and BOD special, Respiration/Dilution, Soil respiration, Biodegradability, Biogas determination |
| BOD range | 0 –4,000 mg/l | 0 –400.000 mg/l |
| Measured value memory | 1 –7 days | 0.5 h – 180 days |
| Pressure mode | — | Pressure p 500 – 1.350 hPa (1500 hPa with OxiTop®-IDS/B measuring head) |
| Sample volume | Fixed | Fixed and definable |

| Composition of complete packages | | | | | |
|--------------------------------------------------------------------------------|---------------------------------------------|------------------------------------------------------------------------|------------------------------------------------------------|--------------------------------------------------------------------------|------------------------------------|
| | | | | | |
| | OxiTop®-i | OxiTop®-IDS | | | |
| Accessories | IS 6/ IS 12 | SET6 / 12 | B6 / B6M / B6M 2.5 | A6 / A12 | AN6 / AN12 |
| Vessels including adapters | Bottle PF 600, amber, 510 ml, rubber sleeve | Bottle PF 600, amber, 510 ml, rubber sleeve | Bottle PF45/500 with adapter AD/SK / vessel MG 1.0/ MG 2.5 | Bottle PF45/1000 with adapter AD/SK / Bottle PF45/250 with adapter AD/SK | Bottle MF45/1000 / Bottle MF45/250 |
| Units | 6 / 12 | 6 / 12 | 6 / 6 / 6 | 6 / 12 | 6 / 12 |
| Measuring heads | OxiTop®-i | OxiTop®-IDS | OxiTop®-IDS | OxiTop®-IDS | OxiTop®-IDS/B |
| Stirring platforms | IS 6/IS 12 | IS 6/IS 12 | — | IS 6-Var/IS 12 | IS 6-Var/IS 12 |
| IDS Multi-parameter instrument including IDSM WLM-M module, software and cable | — | Multi 3620 IDS Set WL/Multi 3630 IDS SETWL, please order separately | | | |
| CO ₂ Absorbent | | | | | |
| Nitrification inhibitor | | | — | | |
| Overflow measuring flasks | 164/432 ml | 164/432 ml | — | — | — |
| Stirring bars | 6/12 | 6/12 | — | 6/12 | 6/12 |
| Stirring bar remover | | | — | | |
| 1 year guaranty for defects of quality on all sensor heads. | | | | | |

| Model | Description | Order No. |
|------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| Benchtop Meters for BOD/Respiration | | |
| inoLab@Oxi 7310 | Professional dissolved oxygen meter, menu controlled, for measurements / documentation according GLP/AQA. For AC and battery operation. Meter with universal power supply, stand and operation manual, software and USB cable. | 1BA300 |
| inoLab@Oxi 7310P | Professional dissolved oxygen meter, menu controlled, for measurements / documentation according GLP/AQA, with integrated thermal printer. Meter with universal power supply, stand and operation manual, software and USB cable. | 1BA300P |
| inoLab@Oxi 7310 SET1 | Professional dissolved oxygen meter, menu controlled, for measurements / documentation according GLP/AQA. Set including galvanic D.O. sensor. Meter with universal power supply, stand and operation manual, galvanical D.O. sensor CellOx@325, abrasive foil, electrolyte, cleaning solution, replacement caps, CD-ROM incl. Software, USB cable. | 1BA301 |
| inoLab@Oxi 7310 SET4 | Professional dissolved oxygen meter, menu controlled, for measurements / documentation according GLP/AQA. Set including galvanic D.O. sensor. Meter with universal power supply, stand and operation manual, self stirring galvanic D.O. sensor StirrOx@G, abrasive foil, electrolyte, cleaning solution, replacement caps, CD-ROM incl. Software, USB cable. | 1BA304 |
| BOD Sensors and Accessories | | |
| StirrOx@G | Self-stirring oxygen sensor for determination in Karlsruhe bottles and Winkler bottles with OxiCal@-ST calibration and storage vessel, 3 exchange membrane heads, cleaning solution, electrolyte and abrasive foil. | 201425 |
| KF 12 | Karlsruhe bottles, incl. 60 mm ground stoppers suitable for BOD determination. Minimum order 10 bottles. | 205700 |
| TS 19 | Funnel set for BOD determination in NS 19 Winkler bottles. | 205710 |
| FD 19/KA | Fixing cover for Karlsruhe bottle KF 12. | 205711 |
| FR 19 | Fixing ring, recommended for measurements in Karlsruhe bottles for CellOx@325 or FDO@925 | 205712 |
| RZ 300 | Attachment stirrer for constant approach flow in combination with a magnetic stirrer. For CellOx@325, FDO@925. | 203824 |
| OxiTop@-i Self Check (effective from February 2019) | | |
| OxiTop@-i IS 6 | <i>New</i> BOD instrument for self-check with 6 measuring units, with direct display of the BOD value at the measuring head, universal power supply 100 to 240 VAC, 50/60 Hz; accessories: see table Composition of complete packages. | 208260 |
| OxiTop@-i IS 12 | <i>New</i> BOD instrument for self-check with 12 measuring units, with direct display of the BOD value at the measuring head, universal power supply 100 to 240 VAC, 50/60 Hz; accessories: see table „Composition of complete packages“. | 208262 |
| OxiTop@-IDS for Self Check (effective from February 2019) | | |
| OxiTop@-IDS Set 6 | <i>New</i> Respirometric measuring system for aerobic measuring, also for BOD self check measurements, 6 measuring units, requires Multi 3620 IDS SETWL or Multi 3630 IDS SETWL, please order separately. Universal power supply 100 to 240 VAC, 50/60 Hz. Accessories: see table „Composition of complete packages“, incl. stirring platform IS6. | 209710 |
| OxiTop@-IDS Set 12 | <i>New</i> Respirometric measuring system for aerobic measuring, also for BOD self check measurements, 12 measuring units, requires Multi 3620 IDS SETWL or Multi 3630 IDS SETWL, please order separately. Universal power supply 100 to 240 VAC, 50/60 Hz. Accessories: see table „Composition of complete packages“, incl. stirring platform IS12. | 209712 |
| OxiTop@-IDS IS 6 | <i>New</i> Extension set for respirometric aerobic OxiTop@-IDS measurement systems (also for BOD self check measurements), stirring platform IS6 for 6 measuring units, universal power supply 100 to 240 VAC, 50/60 Hz; accessories: including 6 sample bottles, 6 OxiTop@-IDS measuring heads, 6 rubber sleeves and 6 magnetic stirring bars. | 209714 |
| OxiTop@-IDS IS 12 | <i>New</i> Extension set for respirometric aerobic OxiTop@-IDS measurement systems (also for BOD self check measurements), stirring platform IS12 for 12 measuring units, universal power supply 100 to 240 VAC, 50/60 Hz; accessories: including 12 sample bottles, 12 OxiTop@-IDS measuring heads, 12 rubber sleeves and 12 magnetic stirrer bars. | 209716 |
| OxiTop@-i Accessories | | |
| OxiTop@-i SET2 | <i>New</i> SET for 2 measuring units, complete with 2 OxiTop@-i measuring heads (1 x blue, 1 x grey), including 2 sample bottles PF600, 2 rubber sleeves GK 600 and 2 magnetic stirring bars RST600. | 208264 |
| OxiTop@-i SET6 | <i>New</i> SET for 6 measuring units, complete with 6 OxiTop@-i measuring heads (3 x blue, 3 x grey), including 6 sample bottles PF600, 6 rubber sleeves GK 600 and 6 magnetic stirring bars RST600. | 208266 |
| OxiTop@-i G | <i>New</i> OxiTop@-i measuring head (grey), menu-operated, with backlit graphic LCD, control LED, operating buttons, for direct BOD readings, volume pre-selection, measurement duration between 1 and 7 days selectable, display of measuring curve, display of current readings. | 208270 |
| OxiTop@-i B | <i>New</i> OxiTop@-i measuring head blue, menu-operated, with backlit graphic LCD, control LED, operating buttons, for direct BOD readings, volume pre-selection, measurement duration between 1 and 7 days selectable, display of measuring curve, display of current readings. | 208272 |
| OxiTop@-i 2 | <i>New</i> 2 OxiTop@-i measuring heads (blue and grey), menu-operated, with backlit graphic LCD, control LED, operating buttons, for direct BOD readings, volume pre-selection, measurement duration between 1 and 7 days selectable, display of measuring curve, display of current readings. | 208274 |

| Model | Description | Order No. |
|------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| OxiTop®-i 6 | <i>New</i> 6 OxiTop®-i measuring heads (3 x blue and 3 x grey), menu-operated, with backlit graphic LCD, control LED, operating buttons, for direct BOD readings, volume pre-selection, measurement duration between 1 and 7 days selectable, display of measuring curve, display of current readings. | 208276 |
| OxiTop®-i 12 | <i>New</i> 12 OxiTop®-i measuring heads (6 x blue and 6 x grey), menu-operated, with backlit graphic LCD, control LED, operating buttons, for direct BOD readings, volume pre-selection, measurement duration between 1 and 7 days selectable, display of measuring curve, display of current readings. | 208278 |
| OxiTop®-IDS Accessories | | |
| IDS WLM-M | Wireless module for connecting to MultiLine®3310/3510/36x0 IDS, inoLab®Multi IDS and MonoLine. Connects up to three sensors at the same time (depends on meter capabilities). Also for operation of OxiTop®-IDS. | 108142 |
| OxiTop®-IDS | <i>New</i> OxiTop®-IDS wireless measuring head with Bluetooth®LE technology for operation with Multi 3620 IDS or Multi 3630 IDS, memory for 360 measuring points, with backlit graphic LCD and control LED, for routine BOD also operable without IDS multi-parameter instrument. | 209700 |
| OxiTop®-IDS/B | <i>New</i> OxiTop®-IDS/B wireless measuring head with Bluetooth®LE technology for operation with Multi 3620 IDS or Multi 3630 IDS, memory for 360 measuring points, with backlit graphic LCD and control LED, enhanced pressure range up to 1.5 bars, H ₂ S resistant pressure sensor for biogas measurement, for routine BOD also operable without IDS multi-parameter instrument. | 209701 |
| OxiTop®-IDS 2 | <i>New</i> Two OxiTop®-IDS wireless measuring heads with Bluetooth®LE technology for operation with Multi 3620 IDS or Multi 3630 IDS, memory for 360 measuring points, with backlit graphic LCD and control LED, for routine BOD also operable without IDS multi-parameter instrument. | 209702 |
| OxiTop®-IDS 6 | <i>New</i> Six OxiTop®-IDS wireless measuring heads with Bluetooth®LE technology for operation with Multi 3620 IDS or Multi 3630 IDS, memory for 360 measuring points, with backlit graphic LCD and control LED, for routine BOD also operable without IDS multiparameter instrument. | 209704 |
| Multi 3620 IDS SETWL | <i>New</i> Professional digital multi meter for portable field measurement, with dual channel input, wireless ready, color graphic display including data logger and USB interfaces, also for controlling OxiTop®-IDS. Instrument with rechargeable batteries, driver software for USB, cable and universal power supply, protection SM Pro, and wireless meter module IDSWLM-M. | 2FD56W |
| Multi 3630 IDS SETWL | <i>New</i> Professional digital multi meter for portable field measurement, with triple channel input, wireless ready, color graphic display incl. data logger and USB interfaces, also for controlling OxiTop®-IDS. Instrument with rechargeable batteries, driver software for USB, cable and universal power supply, protection SM Pro, and wireless meter module IDSWLM-M. | 2FD57W |
| Stirring Platforms | | |
| IS 6 | Stirrer for 6 measuring units, universal power supply 100 to 240 VAC, 50/60 Hz; | 208140 |
| IS 12 | Stirrer for 12 measuring units, universal power supply 100 to 240 VAC, 50/60 Hz; | 208141 |
| IS 6-Var | Stirrer for 6 sample or measuring vessels PF../MF../MG.. with base diameter up to 93 mm, universal power supply 100 to 240 VAC, 50/60 Hz | 208145 |
| OxiTop®-IDS for Determination of Soil Respiration (effective from February 2019) | | |
| OxiTop®-IDS B6 | <i>New</i> Measuring system for determination of soil respiration (aerobic), 6 measuring units, requires Multi 3620 IDS SETWL or Multi 3630 IDS SETWL, please order separately. Accessories: see table „Composition of complete packages“, incl. 6 bottles PF45/500 of 500 ml each and 6 OxiTop®AD/SK adapters. | 209728 |
| OxiTop®-IDS B6M | <i>New</i> Measuring system for determination of soil respiration AT4 (aerobic), 6 measuring units with 2.5 l measuring vessels, requires Multi 3620 IDS SETWL or Multi 3630 IDS SETWL, please order separately. Accessories: see table „Composition of complete packages“, including 6 vessels MG 2.5, each including lid adapter DV/MG, DV/MG/KL clamps, GK 600 rubber sleeve and beaker MB 50, CO ₂ absorbent. | 209730 |
| OxiTop®-IDS B6M-2.5 | <i>New</i> Measuring system for determination of soil respiration according AT4 (aerobic) with IDS multiparameter instrument Multi 3620 IDS and 6 measuring units with 2.5 l measuring vessels. Accessories: see table „Composition of complete packages“, including 6 vessels MG 2.5, each including lid adapter DV/MG, DV/MG/KL clamps, GK 600 rubber sleeve and beaker MB 50, CO ₂ absorbent. | 209732 |
| OxiTop®-IDS for Determination of Biodegradability acc. To OECD (effective from February 2019) | | |
| OxiTop®-IDS Set 6 | <i>New</i> Respirometric measuring system for aerobic measuring, also for BOD self check measurements, 6 measuring units, requires Multi 3620 IDS SETWL or Multi 3630 IDS SETWL, please order separately. Universal power supply 100 to 240 VAC, 50/60 Hz. Accessories: see table „Composition of complete packages“, incl. stirring platform IS6. | 209710 |
| OxiTop®-IDS Set 12 | <i>New</i> Respirometric measuring system for aerobic measuring, also for BOD self check measurements, 12 measuring units, requires Multi 3620 IDS SETWL or Multi 3630 IDS SETWL, please order separately. Universal power supply 100 to 240 VAC, 50/60 Hz. Accessories: see table „Composition of complete packages“, incl. stirring platform IS12. | 209712 |
| OxiTop®-IDS IS 6 | <i>New</i> Extension set for respirometric aerobic OxiTop®-IDS measurement systems (also for BOD self check measurements), stirring platform IS6 for 6 measuring units, universal power supply 100 to 240 VAC, 50/60 Hz; accessories: including 6 sample bottles, 6 OxiTop®-IDS measuring heads, 6 rubber sleeves and 6 magnetic stirring bars. | 209714 |

| Model | Description | Order No. |
|-----------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| OxiTop®-IDS IS 12 | <i>New</i> Extension set for respirometric aerobic OxiTop®-IDS measurement systems (also for BOD self check measurements), stirring platform IS12 for 12 measuring units, universal power supply 100 to 240 VAC, 50/60 Hz; accessories: including 12 sample bottles, 12 OxiTop®-IDS measuring heads, 12 rubber sleeves and 12 magnetic stirrer bars. | 209716 |
| OxiTop®-IDS A6 | <i>New</i> Respirometric measuring system for aerobic operation, 6 measuring units, requires Multi 3620 IDS SETWL or Multi 3630 IDS SETWL, please order separately. Universal power supply 100 to 240 VAC, 50/60 Hz. Accessories: see table „Composition of complete packages“, incl. 6 sample bottles, PF45/1000 (1000 ml), 6 OxiTop®AD/SK adapters and 1 stirring platform IS 6-Var. | 209720 |
| OxiTop®-IDS A12 | <i>New</i> Respirometric measuring system for aerobic operation, 12 measuring units, requires Multi 3620 IDS SETWL or Multi 3630 IDS SETWL, please order separately. Universal power supply 100 to 240 VAC, 50/60 Hz. Accessories: see table „Composition of complete packages“, incl. 12 sample bottles, PF45/250 (250 ml), 12 OxiTop®AD/SK adapters and 1 stirring platform IS 12. | 209722 |
| OxiTop®-IDS for Biogas Determination/Microbiology (effective from February 2019) | | |
| OxiTop®-IDS AN6 | <i>New</i> Respirometric measuring system for anaerobic operation, 6 measuring units, requires Multi 3620 IDS SETWL or Multi 3630 IDS SETWL, please order separately. Universal power supply 100 to 240 VAC, 50/60 Hz. Accessories: see table „Composition of complete packages“, incl. 6 MF45/1000 bottles (1000 ml) with 2 nozzles each and septum seal, 6 OxiTop®AD/SK adapters and 1 stirring platform IS 6-Var. | 209724 |
| OxiTop®-IDS AN12 | <i>New</i> Respirometric measuring system for anaerobic operation with IDS multiparameter instrument Multi 3620 IDS and 12 measuring units, universal power supply 100 to 240 VAC, 50/60 Hz. universal power supply 100 to 240 VAC, 50/60 Hz. Accessories: see table „Composition of complete packages“, incl. 12 MF45/250 bottles (250 ml) with 2 nozzles each and septum seal, 12 OxiTop®AD/SK adapters and 1 stirring platform IS 12. | 209726 |
| Accessories/Spare Parts Respirometry | | |
| OxiTop®AD/SK | Adapter for OxiTop®/ OxiTop®-C measuring heads for all sample and measuring bottles models PF45/..and MF/45, suitable for autoclave (incl. OxiTop®AD adapter, OxiTop®SK screwable adapter, gasket DI 45 for adapter, screwing set VS45 for adapter, 2x GK 600 rubber sleeves). | 209016 |
| MG 1,0 | Measuring vessel 1000 ml for OxiTop®-C. Accessories: 1x DV/MG adapter with rack and gasket, 2x rubber sleeves GK 600, 4x DV/MG/KL clamps and 1x MB 50 measuring beaker (50 ml). | 209019 |
| DV/MG | Cover adapter for measuring vessels MG 1.0 and MG 1.5 complete with holder and seal DV/MG/DI. | 209022 |
| MG 2.5 | Measuring vessel for OxiTop®C with 2.5 l volume, accessories 1x V/MG lid adapter with rack and seal, 2x GK 600 rubber sleeve, 4x DV/MG/KL clamps and beaker MB 50 (50 ml). | 209023 |
| PF45/250 | Sample bottle, contents: 250 ml; Minimum order 3 pieces (packing unit), price per piece. | 209102 |
| PF45/500 | Sample bottle, contents: 500 ml; Minimum order 3 pieces (packing unit), price per piece. | 209103 |
| PF45/1000 | Sample bottle, contents: 1000 ml; Minimum order 3 pieces (packing unit), price per piece. | 209104 |
| MF45/250 | Measuring bottle, contents: 250 ml with 2 nozzles incl. 2x VS25 screwing sets for nozzle, 1x GK 600 rubber sleeve, 2x SEP25 Septum for nozzle 2x STH25 support for nozzle. | 209106 |
| MF45/500 | Measuring bottle, contents: 500 ml with 2 nozzles incl. 2x VS25 screwing sets for nozzle, 1x GK 600 rubber sleeve, 2x SEP25 Septum for nozzles, 2x STH25 support sleeves for nozzle. | 209107 |
| MF45/1000 | Measuring bottle, contents: 1000 ml with 2 nozzles incl. 2x VS25 screwing sets, 2x SEP25 Septum for nozzles, 2x STH25 support sleeves for nozzle. | 209108 |
| MF45/250/E | Spare measuring bottle, contents 250 ml. | 209110 |
| MF45/500/E | Spare measuring bottle, contents 500 ml. | 209111 |
| MF45/1000/E | Spare measuring bottle, cont. 1000 ml. | 209112 |
| NKI 500 | Soda lime as CO ₂ absorbant, with indicator, 500 g. | 209142 |
| VS 45 | Screw connection for adapter OxiTop®AD/SK /OxiTop®AD; minimum order 3 pieces (packing unit), price per piece. | 209500 |
| DI 45 | Seal for OxiTop®AD/SK / OxiTop®AD adapter suitable for autoclave (121°C); minimum order 5 pieces (packing unit), price per piece. | 209501 |
| VS 25 | Screwable connection for nozzle of measuring vessels MF45/..., suitable for autoclave, minimum order 5 pieces (packing unit), price per piece. | 209502 |
| SEP 25 | Septum for nozzle of measuring vessels MF45/..., suitable for autoclave (121°C), minimum order 5 pieces (packing unit), price per piece. | 209503 |
| STH 25 | Supporting sleeve for nozzle of measuring vessels MF45/..., suitable for autoclave (121°C), minimum order 5 pieces (packing unit), price per piece. | 209504 |
| MG 1,0 E | Replacement measuring vessel 1000 ml, minimum order 6 pcs, price per piece. | 209510 |
| MG 2.5 E | Replacement vessel MG 2.5, minimum order 4 units, price per unit. | 209512 |
| DV/MG/KL | Clip for measuring vessel MG 1.0, minimum order 20 pcs, price per piece. | 209520 |
| DV/MG/DI | Seal for measuring vessel MG 1.0 and MG 2.5, minimum order 5 pcs, price per piece. | 209521 |
| MB 50 | Beaker, 50 ml | 108064 |

| Model | Description | Order No. |
|--------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| BOD Thermostat Boxes | | |
| OxiTop® Box | BOD OxiTop® thermostat box, for max. 12 OxiTop® units/20 Karlsruhe bottles for IS/6 / IS/12 stirrer, 230 VAC, 50 Hz. | 208432 |
| OxiTop® Box (115 VAC) | BOD OxiTop® thermostat box, for max. 12 OxiTop® units/20 Karlsruhe bottles for IS/6 / IS/12 stirrer, 115 VAC, 60 Hz. | 208433 |
| Thermostat Cabinets | | |
| AK IS6/TS 606/2i | Replacement Y-cable for operating stirring platforms with TS 606/2-i respectively TS 606-G/2-i effective from December 2013. | 208386 |
| TS 608/2-i | <i>New</i> Incubator cabinet for 2 BOD OxiTop® systems (12 standard bottles each), 230 VAC, 50 Hz, refrigerant R-600a. | 208450 |
| TS 608-G/2-i | <i>New</i> Incubator cabinet for 2 BOD OxiTop® systems (12 standard bottles each), with insulating glass front door, 230 VAC, 50 Hz, refrigerant R-600a. | 208452 |
| TS 608/4-i | <i>New</i> Incubator cabinet for 4 BOD OxiTop® systems (12 standard bottles each), 230 VAC, 50 Hz, refrigerant R-600a. | 208454 |
| TS 608-G/4-i | <i>New</i> Incubator cabinet for 4 BOD OxiTop® systems (12 standard bottles each), with insulating glass front door, 230 VAC, 50 Hz, refrigerant R-600a. | 208456 |
| TS 1008-i | <i>New</i> Incubator cabinet for 4 BOD OxiTop® systems, special applications with large compartment height for 1.5 l flasks/bottles with nozzle, 230 VAC, 50 Hz, refrigerant R-600a. | 208458 |
| Additional Accessories BOD | | |
| ST-OxiTop® | Storage rack for 6 OxiTop® or 6 OxiTop®-C measuring heads. | 209010 |
| Mark-6 | Set of marker rings (6 pcs) and identification shields (1–12) for identification of BOD bottles for OxiTop® systems. | 209013 |
| Dia/600 | 5 blocks of chart paper BOD. | 209020 |
| MK 22/600 | Overflow measuring flask, contents 22.7 ml | 209030 |
| MK 43/600 | Overflow measuring flask, contents 43.5 ml | 209040 |
| MK 97/600 | Overflow measuring flask, contents 97.0 ml | 209050 |
| MK 164/600 | Overflow measuring flask, contents 164 ml | 209060 |
| MK 250/600 | Overflow measuring flask, contents 250 ml | 209070 |
| MK 365/600 | Overflow measuring flask, contents 365 ml | 209080 |
| MK 432/600 | Overflow measuring flask, contents 432 ml | 209090 |
| PF 600 | Sample bottle, amber, contents 510 ml, minimum order 3 pcs, price per piece. | 209100 |
| RST 600 | Stirrer bar, minimum order 3 pcs, price per piece. | 209120 |
| REF 600 | Stirrer bar remover. | 209130 |
| NHP 600 | Sodium hydroxide pellets, 2 bottles of 50 g each. | 209140 |
| GK 600 L | Rubber sleeve L for OxiTop® systems usable with PF 600 and MG 1,0 with DV/MG adapter. | 209172 |
| NTH 600 | Nitrification inhibitor, light protected bottle with content 25 ml. Concentration 5 g/l allylthiourea | 209331 |
| OxiTop® PM | System check of OxiTop® systems with 10 calibration tablets for checking measuring units (long term test of leakage testing and negative pressure). | 209333 |
| OxiTop® PT | Check unit for negative pressure and leakage testing for all OxiTop® measuring heads. | 209334 |
| OxiTop®-C Accessories (old systems) | | |
| OxiTop®-C 2 | 2 OxiTop®-C measuring heads for use with OxiTop®-Controller OC 100/OC 110 with IR interface and data storage for 360 data sets. | 208823 |
| OxiTop®-C | OxiTop®-C replacement measuring head for use with OxiTop®-Controller OC 100/OC 110 with IR interface and data storage for 360 data sets. | 208830 |
| OxiTop®-C/B | OxiTop®-C/B replacement measuring head for use with OxiTop®-Controller OC 100/OC 110 with IR interface and memory for 360 data sets. Corrosion proof model for anaerobic/biogas applications. | 208831 |
| AK 540/B | Interface cable for connection of instrument with 6-pin socket to a PC via RS232 interface. | 902842 |

Distributed By: **Camlab Ltd**
Unit 24, Norman Way Industrial Estate
Over, Cambridge, CB24 5WE, United Kingdom
T: +44 (0) 1954 233 110 E: sales@camlab.co.uk



