

Method and Reagents for CW3000 and CW5000

| | Measuring Range | Chemistry | Reagent Type | Pack Size | Reference | Price | CW3000 | CW5000 |
|------------------|-----------------|--------------------------|----------------------------|-----------|---------------|--------|--------|--------|
| Alkalinity-m | 5-200mg/l | Acid/Indicator | Tablet/Indicator | 100 | TT/51.32.10 | £9.70 | • | • |
| Alkalinity-p | 5-500mg/l | Acid/Indicator | Tablet/Indicator | 100 | TT/51.32.30 | £13.60 | • | • |
| Aluminium | 0.01-0.25mg/l | Eriochrome cyanine R | Pillow/Indicator | Set | CW/53.50.00 | £71.00 | • | • |
| Ammonia | 0.02-1mg/l | Salicylate | Tablet Ammonia No1 | 100 | TT/51.25.80 | £9.00 | • | • |
| | | | Tablet Ammonia No2 | 100 | TT/51.25.90 | £9.00 | • | • |
| Ammonia | 0- 0.5mg/l | Salicylate | Powder/ liquid | 100 | CW/53.55.00 | £83.50 | | • |
| Ammonia | 0-2.5mg/l | Salicylate | Tube | Set | CW/53.56.00 | £81.65 | • | • |
| Ammonia | 0-50mg/l | Salicylate | Tube | Set | CW/53.56.50 | £81.65 | • | • |
| Bromine | 0.05-6.5mg/l | DPD | Tablet DPD No1 | 100 | TT/51.10.60 | £6.50 | • | • |
| | | | Tablet DPD No1 High Ca | 100 | TT/51.57.40 | £7.70 | • | • |
| Bromine | 0.1-4.5mg/l | DPD | Liquid DPD 1 Buffer | 15ml | TT/47.10.10 | £10.90 | | • |
| | | | Liquid DPD 1 Reagent | 15ml | TT/47.10.20 | £10.90 | | • |
| Cadmium | 0.025-0.75mg/l | Cadion | Spectroquant | 25 | TT/42.07.50 | | | • |
| Chloride | 0.5-25mg/l | Silver nitrate/turbidity | Tablet Chloride T1 | 100 | TT/51.59.10 | £10.90 | • | |
| | | | Tablet Chloride T2 | 100 | TT/51.59.20 | £10.90 | • | |
| Chlorine | 0.02-3mg/l | DPD | Tablet DPD No1 | 100 | TT/51.10.60 | £6.50 | • | • |
| | | | Tablet DPD No1 High Ca | 100 | TT/51.57.40 | £7.70 | • | • |
| | | | Tablet DPD No3 | 100 | TT/51.10.80 | £6.70 | • | • |
| Chlorine | 0.05-2mg/l | DPD | Liquid DPD 1 Buffer | 15ml | TT/47.10.10 | £10.90 | • | • |
| | | | Liquid DPD 1 Reagent | 15ml | TT/47.10.20 | £10.90 | • | • |
| | | | Liquid DPD 3 | 15ml | TT/47.10.30 | £10.90 | • | • |
| Chlorine | 0.01-2mg/l | DPD | Pillow DPD free | 100 | CW/53.01.00 | £11.50 | • | • |
| | | | Pillow DPD total | 100 | CW/53.01.20 | £11.50 | • | • |
| Chlorine HR | 5-200mg/l | Potassium Iodide | Tablet Chlorine HR(KI) | 100 | TT/51.30.00 | £12.40 | • | • |
| | | | Tablet Acidifying GP | 100 | TT/51.54.80 | £9.20 | • | • |
| Chlorine Dioxide | 0.05-11mg/l | DPD | Tablet DPD No1 | 100 | TT/51.10.60 | £6.50 | • | |
| | | | Tablet DPD No3 | 100 | TT/51.10.80 | £6.70 | • | |
| | | | Tablet Glycine | 100 | TT/51.21.70 | £8.50 | • | |
| Chlorine Dioxide | 0.04-2.5mg/l | DPD | Tablet DPD No1 | 100 | TT/51.10.60 | £6.50 | | • |
| | | | Tablet DPD No1 High Ca | 100 | TT/51.57.40 | £7.70 | | • |
| Chlorine Dioxide | 0.5-2.5mg/l | DPD | Liquid DPD 1 Buffer | 15ml | TT/47.10.10 | £10.90 | | • |
| | | | Liquid DPD 1 Reagent | 15ml | TT/47.10.20 | £10.90 | | • |
| Chromium | 0.005-2mg/l | 1,5-Diphenyl-carbazone | Pillow Persulphate | 100 | CW/53.73.00 | £28.40 | | • |
| | | | Pillow Chromium Hexavalent | 100 | CW/53.73.10 | £21.20 | | • |
| COD | 0-150mg/l | Dichromate/H2SO4 | COD Vial | 25 | CW/2.42.07.20 | £27.70 | • | • |
| | 0-1500mg/l | Dichromate/H2SO4 | COD Vial | 25 | CW/2.42.07.21 | £27.70 | • | • |
| | 0-15000mg/l | Dichromate/H2SO4 | COD Vial | 25 | CW/2.42.07.22 | £27.70 | • | • |
| Copper | 0.05-5mg/l | Biquinoline | Tablet Copper No1 | 100 | TT/51.35.50 | £27.70 | • | • |
| | | | Tablet Copper No2 | 100 | TT/51.35.60 | £11.90 | • | • |