

Method and Reagents for CW3000 and CW5000

	Measuring Range	Chemistry	Reagent Type	Pack Size	Reference	Price	CW3000	CW5000
Copper	0.05-5mg/l	Bicinchoninate	Pillow Copper 1	100	CW/53.03.00	£23.60	•	
Cyanide	0.005-0.5mg/l	Pyridine-barbituric acid	Powder/Liquid	Set	TT/2.41.88.75		•	•
Cyanuric Acid	2-160mg/l	Melamine	Tablet Cyanuric Acid	100	TT/51.13.20	£12.30	•	
DEHA	20-500ug/l	PPST*	DEHA Solution	100ml	TT/46.11.81	£7.40	•	•
			Tablet DEHA	100	TT/51.32.20	£42.30	•	•
DEHA	20-500ug/l	PPST*	Powder/liquid	100	CW/53.60.00	£50.00	•	•
Fluoride	0.02-1.5mg/l	SPADNS	Liquid Reagent SPADNS	100ml	TT/46.74.81	£25.50	•	•
Formaldehyde	0.1-5mg/l	H ₂ SO ₄ /Chromotropic Acid	Tube	25	TT/42.07.52	£83.40		•
Hardness, Calcium	50-900mg/l	Murexide	Tablet Calcheck	100	TT/51.56.50	£10.40	•	
Hardness, Total	2-50mg/l	Metallphthalein	Tablet Hardcheck	100	TT/51.56.60	£10.40	•	•
Hazen	0-500Pt/Co units	Direct Reading					•	
Hydrogen Peroxide	0.01-1.5mg/l	DPD/Catalyst	Tablet Hydrogen Peroxide	100	TT/51.23.80	£12.30	•	•
Iodine	0.05-3.6mg/l	DPD	Tablet DPD No1	100	TT/51.10.60	£6.50	•	•
Iron Soluble	0.01-1mg/l	PPST	Tablet Iron LR	100	TT/51.53.70	£21.00	•	•
			Tablet Iron LR	100	TT/51.54.20	£21.00	•	•
Iron Fe2+	0.01-1mg/l	PPST	Tablet Iron LR	100	TT/51.53.70	£21.00	•	•
			Tablet Iron LR	100	TT/51.54.20	£21.00	•	•
Iron Soluble	0-3mg/l	1,10-Phenanthroline	Pillow Vario Ferro	100	CW/53.05.60	£12.10	•	•
Iron TPTZ	0.02-1.8mg/l	TPTZ	Pillow Vario TPTZ	100	CW/53.05.50	£16.60	•	•
Lead	0.1-5mg/l	4-(2-Pyridylazo)-resorcin	Tube	25	TT/42.07.54	£83.40		•
Manganese	0.05-4mg/l	Formaloxin	Tablet Manganese 1	100	TT/51.60.80	£28.30	•	•
			Tablet Manganese 2	100	TT/51.60.90	£12.10	•	•
Manganese LR	0.01-0.7mg/l	PAN	Pillow/Liquid	Set	CW/53.50.90	£46.50	•	•
Molybdate	1-50mg/l	Thioglycolate	Tablet Molybdate No1	100	TT/51.30.60	£9.60	•	•
			Tablet Molybdate No2	100	TT/51.30.70	£9.60	•	•
Nickel	0.02-7mg/l	Dimethylglyoxime	Tube	25	TT/2.41.90.33	£36.20		•
Nitrate	0-30mg/l	Chromotropic Acid	Powder/Tube	Set	CW/53.55.80	£84.50	•	•
Nitrite	0.01-0.5mg/l	NED	Tablet Nitrite LR	100	TT/51.23.10	£9.00	•	•
Nitrite	0-3mg/l	Sulfanilic Acid/Naphthyl-amine	Tube	24	TT/2.41.90.18			•
Nitrogen Total	0.5-25mg/l	Persulphate Digestion	Pillow/Tube	Set	CW/53.55.50	£104.50	•	•
Nitrogen Total	5-150mg/l	Persulphate Digestion	Pillow/Tube	Set	CW/53.55.60	£109.25	•	•
Oxygen Active	0.1-10mg/l	DPD	Tablet DPD No4	100	TT/51.12.20	£8.40	•	
Ozone	0.02-1mg/l	DPD/Glycine	Tablet DPD No1	100	TT/51.10.60	£6.50	•	•
			Tablet DPD No3	100	TT/51.10.80	£6.80	•	•
			Tablet Glycine	100	TT/51.21.70	£8.50	•	•
PHMB	2-60mg/l	Buffer/Indicator	Tablet PHMB	100	TT/51.61.00	£18.20	•	