

DR/800 Parameters and Reagent Sets

| Parameter | Range | No. of tests | Reference | Price | DR/820 | DR/850 | DR/890 |
|--|-------------|--------------|-------------|---------|--------|--------|--------|
| Nitrite HR | 0-150mg/l | 100 | HH/21075-69 | £27.30 | | | • |
| Nitrite LR(Accuvac) | 0-0.35mg/l | 25 | HH/25120-25 | £27.80 | • | • | • |
| Nitrite TNT | 0-0.50mg/l | 50 | HH/26083-45 | £61.30 | • | • | • |
| Nitrogen ammonia | 0-0.50mg/l | 100 | HH/26680-00 | £130.20 | | • | • |
| Nitrogen ammonia LR TNT | 0-2.50mg/l | 50 | HH/26045-45 | £90.20 | | • | • |
| Nitrogen ammonia HR TNT | 0-50mg/l | 50 | HH/26069-45 | £69.60 | | • | • |
| Nitrogen Monochloramine and free ammonia | 0-0.50mg/l | 50 | HH/28797-04 | £62.37 | | • | • |
| Nitrogen TKN | 0-150mg/l | 100 | HH/24593-00 | £58.80 | | | • |
| Nitrogen total | 0-25mg/l | 50 | HH/26722-45 | £115.50 | | | • |
| Nitrogen total HR TNT | 0-150mg/l | 50 | HH/27141-00 | £120.80 | | | • |
| Nitrogen total Inorganic | 0-25mg/l | 50 | HH/26049-45 | £77.40 | | • | • |
| Oxygen DO HR(Accuvac) | 0-15mg/l | 25 | HH/25150-25 | £31.00 | • | • | • |
| Oxygen DO LR(Accuvac) | 0-1000ug/l | 25 | HH/25010-25 | £31.00 | | • | • |
| COD LR | 0-150mg/l | 25 | HH/21258-51 | £31.00 | | | • |
| COD HR | 0-1500mg/l | 25 | HH/21259-51 | £31.00 | | • | • |
| COD HR+ | 0-15000mg/l | 25 | HH/24159-51 | £31.00 | | • | • |
| COD Manganese III | 20-1000mg/l | 25 | HH/26234-51 | £31.00 | • | • | • |
| Ozone LR(Accuvac) | 0-0.25mg/l | 25 | HH/25160-25 | £27.80 | | • | • |
| Ozone MR(Accuvac) | 0-0.75mg/l | 25 | HH/25170-25 | £27.80 | | • | • |
| Ozone HR(Accuvac) | 0-1.5mg/l | 25 | HH/25180-25 | £27.80 | | • | • |
| pH | 6.5-8.5 | 50 | HH/26575-12 | £10.50 | • | • | • |
| Phosphonate | 0-125mg/l | 100 | HH/24297-00 | £92.30 | | • | • |
| Phosphorus(Molybdovanadate) | 0-45mg/l | 100 | HH/20760-32 | £15.20 | | | • |
| Phosphorus (Molybdovanadate-Accuvac) | 0-45mg/l | 25 | HH/25250-25 | £23.20 | | | • |
| Phosphorus(Phosver 3) | 0-2.5mg/l | 100 | HH/21060-69 | £30.50 | | • | • |
| Phosphorus(Phosver 3 - Accuvac) | 0-2.5mg/l | 25 | HH/25080-25 | £27.80 | | • | • |
| Phosphorus(Phosver 3 TNT) | 0-5.0mg/l | 50 | HH/27425-45 | £53.30 | | • | • |
| Phosphorus Total (Phosver 3 TNT) | 0-3.5mg/l | 50 | HH/27426-45 | £93.30 | | • | • |
| Phosphorus Acid Hydrolyzable (Phosver 3 TNT) | 0-5.0mg/l | 50 | HH/27427-45 | ££72.10 | | • | • |
| Phosphorus(Amino acid) | 0-30mg/l | 100 | HH/22441-00 | £23.00 | • | • | • |
| Silica UHR | 0-200mg/l | 100 | HH/22443-00 | £88.60 | | | • |
| Silica HR | 0-75.0mg/l | 100 | HH/24296-00 | £56.10 | | | • |
| Silica LR | 0-1.6mg/l | 100 | HH/24593-00 | £58.80 | | | • |
| Sulphate | 0-70mg/l | 100 | HH/21067-69 | £20.90 | • | • | • |
| Sulphate(Accuvac) | 0-70mg/l | 25 | HH/25090-25 | £27.80 | • | • | • |
| Sulphide | 0-0.7mg/l | 100 | HH/22445-00 | £56.60 | | • | • |
| Suspended solids | 0-750mg/l | - | - | | | • | • |
| Triazole | 0-15mg/l | 100 | HH/21412-99 | £35.70 | | | • |
| Tannin and Lignin | 0-9.0mg/l | 100 | HH/22446-00 | £5.00 | | • | • |
| ToxiTrak Toxicity | %inhibition | 100 | HH/25972-00 | £281.40 | | • | • |
| Turbidity | 1-1000FAU | - | - | | • | • | • |
| TOC HR | 100-700mg/l | 50 | HH/27604-45 | £171.20 | | • | • |
| TOC MR | 15-150mg/l | 50 | HH/28159-45 | £159.00 | | • | • |
| TOC LR | 0-20mg/l | 50 | HH/27603-45 | £171.20 | | • | • |
| Volatile acid (Esterification) | 0-2800mg/l | 90 | HH/22447-00 | £101.30 | • | • | • |
| Zinc | 0-3.0mg/l | 100 | HH/24293-00 | £45.40 | | • | • |

For details on Camlabs own range of powder pillows, TNT and COD vials please call our technical support on 01954 233120 to discuss your requirements.