

# Selection of Digital titrator cartridges and sets available

HACH REAGENTS					
Parameter	Method	Set	Price	Cartridge	Price
Acidity 10-160mg/l	Methyl orange 0.1600 NaOH	HH/22728-00	£58.50	HH/14377-01	£15.30
Acidity 10-160mg/l	Phenolphthalein 0.1600 NaOH	HH/22728-00	£58.50	HH/14377-01	£15.30
Acidity 100-4000mg/l	Methyl orange 1.600 NaOH	HH/22728-00	£58.50	HH/14379-01	£15.30
Alkalinity 10-160mg/l	Phenolphthalein 0.160 H <sub>2</sub> SO <sub>4</sub>	HH/22719-00	£52.30	HH/14388-01	£15.30
Alkalinity 100-4000mg/l	Phenolphthalein 1.60 H <sub>2</sub> SO <sub>4</sub>	HH/22719-00	£52.30	HH/14389-01	£15.30
Alkalinity T 10-160mg/l	Brom Cresol/Methyl Red 0.1600 H <sub>2</sub> SO <sub>4</sub>	HH/22719-00	£52.30	HH/14388-01	£15.30
Alkalinity T 100-4000mg/l	Brom Cresol/Methyl Red 1.600 H <sub>2</sub> SO <sub>4</sub>	HH/22719-00	£52.30	HH/14389-01	£15.30
Carbon dioxide 10-100mg/l CO <sub>2</sub>	Phenolphthalein 0.3636 NaOH	HH/22727-00	£40.60	HH/14378-01	£15.30
Carbon dioxide 100-1000mg/l CO <sub>2</sub>	Phenolphthalein 3.636 NaOH	HH/22727-00	£40.60	HH/14380-01	£15.30
Chloride 10-160mg/l	Mercuric nitrate 0.22556 HgN <sub>2</sub> O <sub>6</sub>	HH/22726-00	£61.30	HH/14393-01	£15.30
Chloride 10-100mg/l	Silver nitrate 0.2256 AgNO <sub>3</sub>	HH/22880-00	£74.30	HH/14396-01	£17.90
Chloride 100-1000mg/l	Silver nitrate 1.128 AgNO <sub>3</sub>	HH/22880-00	£74.30	HH/14397-01	£17.90
Chloride 100-8000mg/l	Mercuric nitrate 2.2556 HgN <sub>2</sub> O <sub>6</sub>	HH/22726-00	£61.30	HH/00921-01	£15.30
Chlorine total 20-2000mg/l	Iodometric 0.113 Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub>	HH/22725-00	£70.50	HH/22673-01	£14.90
Chlorine total 2000-70000mg/l	Iodometric 2.00 Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub>	HH/24448-00	£94.40	HH/14401-01	£14.90
Hardness calcium 10-160mg/l CaCO <sub>3</sub>	EDTA (CalVer) 0.080 EDTA	HH/24472-00	£45.40	HH/14364-01	£15.30
Hardness calcium 100-4000mg/l CaCO <sub>3</sub>	EDTA (CalVer) 0.80 EDTA	HH/24475-00	£49.10	HH/14399-01	£15.30
Hardness total 10-160mg/l CaCO <sub>3</sub>	EDTA (ManVer) 0.080 EDTA	HH/24480-00	£35.20	HH/14364-01	£15.30
Hardness total 100-4000mg/l CaCO <sub>3</sub>	EDTA (ManVer) 0.80 EDTA	HH/24481-00	£35.20	HH/14399-01	£15.30
Hypochlorite 50-150g/l	Iodometric 2026 Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub>	HH/26870-00	£92.80	HH/26869-01	£12.60
Iron 10-100mg/l Fe <sup>3+</sup>	EDTA 0.0716 EDTA	HH/24492-00	£92.80	HH/20817-01	£15.30
Iron 100-1000mg/l Fe <sup>3+</sup>	EDTA 0.716 EDTA	HH/24493-00	£44.80	HH/20818-01	£15.30
Oxygen dissolved 1-10mg/l DO	Mod. Winkler 0.200 NaS <sub>2</sub> O <sub>3</sub>	HH/22722-00	£72.60	HH/22675-01	£15.30
Sulphite 0-800mg/l	Iodate-iodide 0.3998 KI	HH/22723-00	£49.70	HH/14961-01	£15.30