

SPECIFICATIONS

Ranges	-2.00 to +19.99pH, +/-2000mV and -10 to + 110°C
Resolution	0.1/0.01/0.001pH (selectable)
Accuracy	+/-0.002pH
Power	4 x AA type battery
Buffer recognition	Auto 4.01, 6.86, 7.00 and 10.01
Rating	IP67 waterproof submersible to 1 metre

Senslon 1 Portable pH Meter complete with Platinum electrode and instructions.

Reference: HH/51700-10 £414.80

Senslon 1 Portable pH Meter complete with gel-filled electrode and instructions

Reference: HH/51700-11 £365.40

Electrode stand

Reference: HH/45300-00 £96.00

Replacement Batteries, pack of 4

Reference: HH/19380-04 £10.20

Electrode washer

Reference: HH/27047-00 £95.10

LIS Chamber for low ionic strength samples

Reference: HH/51899-00 £88.00

Carrying case (rigid)

Reference: HH/51903-00 £50.70

For full pH electrode listing
see pages 193-198

Hach Senslon 2 Advanced Portable pH/ISE Meter

Features all the attributes of the basic model with the added benefit of being able to download data via the innovative docking station. The IP67 rated casing, back lit display, quick calibration and response and long lasting batteries make it ideal for field work. Will also accept suitable ISE electrodes.

SPECIFICATIONS

Ranges	-2.00 to +19.99pH, -10 to +110°C, 0.000 to 19900 concentration, +/- 2000mv
Resolution	0.1/0.01/0.001pH (selectable)
Accuracy	+/-0.002pH
Power	4 x AA batteries or 220V through docking station
Output	One way RS232 via docking station
Buffer recognition	Auto for 4.01, 6.86, 7.00 and 10.01
Rating	Waterproof IP67, submersible to 1 metre

Senslon 2 Portable pH/ISE meter complete with Platinum electrode and instructions

Reference: HH/51725-10 £630.00



Hach Senslon 2 Advanced
Portable pH/ISE Meter