

# pH



Schott Handylab pH Meter  
Model pH11/pH12

**For full pH electrode listing  
see pages 193-198**



Hanna pH Checker

Senslon 2 Portable pH/ISE meter complete with gel filled electrode and instructions

**Reference:** HH/51725-11 £540.00

Hach Link Software for computer interface

**Reference:** HH/49665-00 £293.00

Docking station

**Reference:** HH/51875-02 £152.30

## Schott Handylab pH Meter Model pH11/pH12

Battery operated meter suitable for the measurement of pH and temperature. Available in two models with the pH12 having the added advantage of data storage and RS232 interface.

### Features:

Range -2.00 to +19.99pH

-5 to +105°C

5000 hours battery life

Low battery indicator

Data storage of 500 records (pH12 model only )

RS232 interface (pH12 model only )

Handylab pH11 portable pH Meter complete with BlueLine23 electrode and instructions

**Reference:** SC/28.520.2917 £390.00

Handylab pH11 portable pH Meter complete with BlueLine24 electrode and instructions

**Reference:** SC/28.520.2982 £430.00

Handylab pH11 portable pH Meter complete with BlueLine14 electrode and instructions

**Reference:** SC/28.520.2999 £450.00

Handylab pH12 portable pH Meter complete with BlueLine24 electrode and instructions

**Reference:** SC/28.520.3054 £530.00

Handylab pH12 portable pH Meter complete with BlueLine14 electrode and instructions

**Reference:** SC/28.520.3062 £540.00

Power supply unit for Handylab12

**Reference:** SC/28.520.5047 £33.00

Connecting cable for PC, for Handylab12

**Reference:** SC/28.520.4959 £142.00

## Hanna pH Checker

The pH checker is a low cost unit with a small slim line electrode allowing measurement in very small containers.

### Features:

- Range 0.00 to 14.00pH

- Resolution 0.01pH

- Accuracy +/-0.2pH

- 2 point calibration

- 3000 hour battery life

- Large LCD display

**Reference:** HI/6602-3650 £27.60

Replacement electrode

**Reference:** HI/6602-3651 £18.90