

pH

For full pH electrode listing
see pages 193-198



Lovibond pH Comparator Kits

Electrode Care and Maintenance Kit

Designed to calibrate, clean and extend the life of your pH electrode. The kit contains pH buffers in twin neck bottles (500ml of each of pH4.00, 7.00 and 10.00 @20°C), 500ml electrode storage solution, electrode cleaning solutions (100ml of each of Biological, Organic and Inorganic), filling solution (100ml of each of 3M KCl/AgCl and 4M KCl), 2 pipettes, 100ml regeneration solution, instruction card and GLP log book.

Reference: RN/RCMK £97.93

Filling Solutions

3mol/litre Potassium chloride (KCl) free from silver ion 100ml

Reference: RN/EFS3 £14.59

3mol/litre Potassium chloride (KCl) free from silver ion 500ml

Reference: RN/EFS35 £36.48

Saturated Potassium Chloride (KCl) free from silver ion 100ml

Reference: RN/LKCL £14.59

Cleaning Solutions

Electrode cleaning solution (Pepsin/hydrochloric acid) for removal of protein 500ml

Reference: RN/ECS £19.45

Inorganic electrode cleaning solution 500ml

Reference: RN/IECS5 £26.75

Organic cleaning solution 500ml

Reference: RN/OECS £21.89

Storage Solutions

pH electrode storage solution 500ml

Reference: RN/ESS5 £21.89

pH electrode storage solution 1000ml

Reference: RN/ESS01 £36.98

Lovibond pH Comparator Kits

A range of easy to use comparator kits complete with reagents and instructions.

Value pH	Volume	Reference	Price
6.5 - 8.4	Phenol red	TT/14.71.00	£29.50
5.2 - 6.8	Bromocresol purple	TT/14.71.10	£29.50
6.0 - 7.6	Bromothymol blue	TT/14.71.20	£29.50
4.0 - 10.0	Universal	TT/14.71.30	£29.50