

Whatman pH Indicator Strips

Whatman Type CF strips have multiple segments of different dye-impregnated paper mounted on each plastic strip. The dyes are chemically bonded so they will not leach out into the solution. Supplied in packs of 100 strips 6 x 80mm.

Whatman CF indicator strips pH 0-14.

Reference: WN/2613-991 £9.00

Whatman CF indicator strips pH 4.5-10

Reference: WN/2614-991 £9.00

Whatman Type CS strips have a printed colour comparison chart on each strip and a central indicator segment which changes colour. This middle is then compared against the chart making it particularly useful for coloured samples. Supplied in packs of 200 strips 11 x 100mm.

Whatman CS indicator strips pH 1-12.

Reference: WN/2612-990 £17.30

Whatman CS indicator strips pH 1.8-3.8.

Reference: WN/2626-990 £14.40

Whatman CS indicator strips pH 3.8-5.5.

Reference: WN/2627-990 £14.40

Whatman CS indicator strips pH 5.2-6.8.

Reference: WN/2628-990 £14.40

Whatman CS indicator strips pH 6.0-8.1.

Reference: WN/2629-990 £14.40

Whatman CS indicator strips pH 8.0-9.7.

Reference: WN/2630-990 £14.40

Whatman pH Indicator Books

Wide range pH indicator book supplied in cartons of 10 packs of 10 books of 20 leaves. The inside cover of each book has a colour comparison chart.

Indicator book pH 1-11

Reference: WN/2600-500 £48.30

Indicator book pH 6.8-8.3.

Reference: WN/2638-500 £40.00



Whatman pH Indicator Strips



Whatman pH Indicator Books