

# Phosphonate/Phosphate

## PHOSPHONATE

### Hach Comparator Kit Model PN-10

Colour disc kit complete with reagents for 100 tests, UV lamp, carrying case and instructions.

**Measuring ranges:** 0 – 5 and 0 -250 mg/l in steps of 1 as  $\text{PO}_4^{3-}$

**Chemistry:** UV photolysis

**Reference:** HH/21133-02 £640.00

Also available without the UV lamp.  
**Reference:** HH/21133-01 £357.00



Hach Comparator Kit

### Hach Pocket Colorimeter

Single parameter handheld colorimeter. Fixed wavelength and three key easy to operate design.

**Measuring ranges:** 0.1 – 2.5mg/l in steps of 0.1 as  $\text{PO}_4^{3-}$   
1 -125mg/l in steps of 0.1 as  $\text{PO}_4^{3-}$

**Chemistry:** UV photolysis

**Features:**

- Silicon photodiode
- LED light source
- IP67 casing
- 2000 test battery life

Supplied with reagents for 100 tests, carrying case, accessories and instructions

**Reference:** HH/58700-08 £540.00



Hach Pocket Colorimeter

## PHOSPHORUS/PHOSPHATE

### Camlab CW2170 Colorimeter

Single parameter handheld colorimeter. Microprocessor controlled for accurate and reproducible results for low range phosphate.

**Measuring range:** 0.05 –4.0mg/l as  $\text{PO}_4$

**Wavelength:** 660nm

**Features:**

- Automatic zero calibration
- LED light source with long term stability
- Large digital display
- Auto switch off

**Reference:** CW/20.79.00 £270.00



Camlab CW2000 Colorimeter