

# Silica/Silver



## Lovibond Comparator Kit

Checkit comparator kit complete with colour disc, reagents and carrying case. For the measurement of low level silica.

**Measuring range:** 0.25 – 4mg/l as SiO<sub>2</sub>

**Reference:** TT/14.73.50 £34.90

## Lovibond Comparator Kit

Checkit comparator kit complete with colour disc, reagents and carrying case. For the measurement of very low level silica.

**Measuring range:** 0 - 1mg/l as SiO<sub>2</sub>

**Reference:** TT/14.73.60 £39.60

## Hach Comparator Kit Model SI-7

Colour disc kit with reagents for 100 tests, carrying case and instructions.

**Measuring range:** 0 – 1mg/l as SiO<sub>2</sub> in steps of 0.02

**Chemistry:** Fast heteropoly blue, amino acid

**Reference:** HH/22550-00 £120.80



## Hach Comparator Kit Model SI-5

Colour disc kit with reagents for 100 tests, carrying case and instructions.

**Measuring ranges:** 0 -40 and 0 - 800mg/l as SiO<sub>2</sub> in steps of 1 and 20 respectively.

**Chemistry:** Fast heteropoly blue, amino acid

**Reference:** HH/14554-00 £112.40

## SILVER

### Hach Rapid Silver test kit

Visual comparison kit measuring in the range 0 – 50ppb as Ag.

Supplied with reagents for 100 tests

**Reference:** HH/26745-00 £182.70

### Silver Ion Selective Electrode

Suitable for any pH/ISE meter, BNC adapters available if required

**For ISE electrodes  
see page 198**