

Tin/Total Organic Carbon (TOC)/Turbidity



TIN

Quantofix Test Papers for Tin

Easy to use test papers supplied in boxes of 100.

Graduations: 0-10-25-50-100-250-500 mg/l Sn²⁺

Reference: MN/913.09

£24.60

TOTAL ORGANIC CARBON (TOC)

For systems to measure TOC please see mutli parameter instruments section.

TOC Standard Solutions (500ml bottles)

Value ppm	Reference	Price
5	RN/TOC5	£41.97
10	RN/TOC10	£41.97
50	RN/TOC50	£41.97
500	RN/TOC500	£41.97
1000	RN/TOC1M	£41.97
2000	RN/TOC2M	£41.97

TURBIDITY

Camlab CW8100 Portable Turbidity Meter

The compact CW8100 portable turbidity meter features an infrared light source and is designed to give fast, precise on-site testing. The unit can be used for clear or coloured solutions and measures the scattered light at an angle of 90° as stipulated in EN ISO7027. The wide measuring range makes the meter suitable for many applications including drinking water, waste water, food and drink and petrochemicals.



SPECIFICATIONS

Measuring range	0.2 to 200 NTU in four ranges
Resolution	0.1NTU in lower two ranges
Measurement cycle	Approximately 9 seconds
Display	LCD
Power supply	9V battery
Ambient conditions	0 - 40°C
Dimensions	190x110x55mm(lxdxh)

Camlab CW8100 Portable Turbidimeter complete with four turbidity standards, battery, samples cells and carrying case.

Reference: CW/20.60.20

£395.00

Set of secondary standards 1, 10, 100 and 1000NTU

Reference: TT/19.36.00

£60.30

Set of 12 sample cells

Reference: TT/19.76.65

£66.30