

# Conductivity



## CONDUCTIVITY

### Camlab Water Model CW6220 Portable Conductivity Meter

The CW6220 microprocessor controlled hand held conductivity meter is designed to meet the day to day demands of reliable on-site and laboratory measurements. The water tight housing rated to IP67 makes it suitable for use in the harshest environments and this is enhanced by the protective housing and built in electrode holder. The easy to use keypad and clear display makes it suitable for field or indoor operation. A clip/stand is incorporated on the rear of the meter to allow for wall hanging or bench use.

#### Features:

- Conductivity, TDS, salinity and temperature.
- Automatic temperature compensation.
- Internal memory for 20 data sets.
- Low battery and battery change indicator.
- Shock-absorbing rubber protective casing.
- Waterproof to IP67.

## SPECIFICATIONS

Conductivity	0.00 to 200.0ms/cm in four ranges
TDS	0 - 1999mg/l
Salinity	0.0 to 70.0g/kg
Temperature	-5.0 to +100.0°C
Power supply	4 x 1.5V battery type AA
Auto-off function	0 - 120 minutes

CW6220 Portable Conductivity Meter complete with measuring cell, accessories, carrying case and instruction manual

**Reference:** CW/72.22.20 £418.00

Spare measuring cell.

**Reference:** TT/72.22.25 £210.00

Replacement carrying case.

**Reference:** TT/72.50.20 £24.40