

Conductivity/Copper

Conductivity Standard Solutions

All supplied in 500ml bottles and calibrated at 25°C

84 uS/cm	
Reference: RN/CSKC84	£18.24
147 uS/cm	
Reference: RN/CSKCS	£18.24
1413 uS/cm	
Reference: RN/CSKCL	£18.24
12880 uS/cm	
Reference: RN/CSKC12880	£18.24
5 uS/cm	
Reference: RN/CSK5	£48.63
10 uS/cm	
Reference: RN/CSKC10	£48.63
50 uS/cm	
Reference: RN/CSKC50	£48.63
100 uS/cm	
Reference: RN/CSKC100	£48.63
500 uS/cm	
Reference: RN/CSKC500	£48.63
1000 uS/cm	
Reference: RN/CSKC1000	£48.63
5000 uS/cm	
Reference: RN/CSK5M	£48.63
10000 uS/cm	
Reference: RN/CSKS10M	£48.63
50000uS/cm	
Reference: RN/CSKS50M	£48.63
100000uS/cm	
Reference: RN/CSKC100M	£48.63
500000uS/cm	
Reference: RN/CSKC500M	£48.63

For other values please contact our **Technical Sales** on **01954 233110**.



Conductivity Standard Solutions

COPPER

Camlab CW2060 Colorimeter

Single parameter handheld colorimeter. Microprocessor controlled for accurate and reproducible results.

Measuring range: 0 – 5mg/l Cu in steps of 0.5

Wavelength: 528nm

Chemistry: Copper No1 and 2

Features:

- Automatic zero calibration
- LED light source with long term stability
- Large digital display
- Auto switch off

Reference: CW/20.72.00 £270.00



Camlab CW2000 Colorimeter