

Cyanide/Cyanuric Acid/DEHA

CYANIDE

Hach Colour Comparator Kit Model CYN-3

Colour disc kit with reagents for 100 tests, instructions and carrying case.

Measuring range: 0 – 0.2 mg/l, CN⁻

Chemistry: Pyridine-Pyrazolone

Reference: HH/02010-02 £165.90

Quantofix Test Papers for Cyanide

Easy to use test papers supplied in boxes of 100.

Graduations: 0-1-3-10-30 mg/l, CN⁻

Reference: MN/913.18 £34.00

Ammonia Ion Selective Electrode

Suitable for any pH/ISE meter,
BNC adapters available if required

For ISE electrodes see page 198



Hach Colour Comparator Kit



Lovibond Mini Kit

CYANURIC ACID

Lovibond Mini Kit

Easy to use turbidity kit

Range 20 – 200mg/l

Reference: TT/41.42.20 £18.40


Hach Colour Comparator Kit Model CY-3

Colour disc kit with reagents for 100 tests, instructions and carrying case.

Measuring range: 20 – 100mg/l.

Chemistry: Turbidity

Reference: HH/01851-02 £63.70



Hach Colour Comparator Kit

DEHA

Camlab CW2070 Colorimeter

Single parameter handheld colorimeter. Microprocessor controlled for accurate and reproducible results.

Measuring range: 0.2 – 0.5mg/l

Wavelength: 528nm

Features:

- Automatic zero calibration
- LED light source with long term stability
- Large digital display
- Auto switch off

Reference: CW/20.97.00 £270.00



Camlab CW2070 Colorimeter