


DEHA/Density



Lovibond Comparator Kit

Lovibond Comparator Kit

Checkit comparator kit complete with colour disc, reagents and carrying case.

Measuring range: 0 – 0.5 mg/l

Chemistry: DEHA reagents

Reference: TT/14.74.30

£34.90



Densito 30PX Portable Density/Specific Gravity Meter

DENSITY

Densito 30PX Portable Density/Specific Gravity Meter

The quick and easy to use Densito 30PX utilises the oscillating tube method to measure the density of liquids. Simply immerse the sample tube into the sample, pull the trigger and take the reading.

Features:

- Small sample size required
- Fast results
- Backlit display which can be activated or deactivated
- Built in time and date
- Instrument identification.
- 1100 results memory

Infrared interface for downloading to PC or printer

SPECIFICATIONS

Range	0-2g/cm ³
Resolution	0.0001g/cm ³
Units	brix%, API with enhanced tables, sulphuric acid %, alcohol-w/w%, alcohol-v/v%, degrees Baume, degrees Plato, and user defined unit.
Accuracy	+/- 0.0001g/cm ³
Temperature Range	+5 to +35°C
Resolution	0.1°C

Densito 30PX portable density meter complete with carry case, sample bottles and CD-ROM

Reference: MT/51324450

£1449.00

Accessories:

Infrared Adapter

Reference: MT/51325006

£159.00

Density water standard sachets 10x6ml

Reference: MT/51325005

£21.00

Syringe adapters

Reference: MT/51324402

£19.00